

INDIVIDUAL KIT™

Designed with all the Tools for the Warfighter/Operator to Mitigate Extremity Bleeding, Tension Pneumothorax and Airway Obstruction.

The Individual Kit from North American Rescue® has been designed for the Warfighter/Operator and packed to provide all of the tools to mitigate the three leading causes of preventable combat death (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).

Our modular system can be attached to the operator using the split-tale connector, PALS/MOLLE style webbing or by fastening to the web/ gun belt using the quick weave method. Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration.

Attachment with split-tail connector and Quick-Grip leg strap or PALS/MOLLE style webbing.

Contents:



Available Colors:



INDIVIDUAL KIT (CCRK)

COLOR	ITEM#	NSN#
BLK	80-0003	6530-01-529-1134
ODG	80-0004	----
COY	80-0005	----

R_X ONLY



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

Features:

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE Style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Available in NAR®'s 11 Point Berry Amendment Nylon Standard

Specifications:

- Case Closed: H 8.5 in. x W 8.75 in. x D 4.5 in.
- Weight: 1 lb 12 oz

Kit Contents:

- 1 x Bag (CCRK®-Individual)
- 4 x Black Talon® Nitrile Trauma Gloves (2 pr.)
- 1 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1 x HyFin® Chest Seal
- 1 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 1 x C-A-T® (Combat Application Tourniquet®)
- 1 x ETD® 6 in. Emergency Trauma Dressing
- 2 x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1 x Polycarbonate Eye Shield
- 1 x Combat Casualty Reference Card
- 1 x Triage Card (NAR T2 Tag®)



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**