

SQUAD KIT

Combat Casualty Response Kits®

CCRK®

REV081623

SQUAD KIT CONTENTS:



Contoured profile, large central compartment, multiple attachment points & front trauma-shear storage

Enhance the Survivability of the Assault Squad

Developed with the end-user requirements and skill set in mind, our innovative Combat Casualty Response Kit® - Squad Kit™ capitalizes on durability, versatility, mobility, and efficiency.

The kit's medical supplies are carefully selected to address the top three causes of death from penetrating trauma as identified by the CoTCCC (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).

Multiple-load carriage configuration capabilities include fanny pack, over the shoulder, attachments for 3-day assault packs and MOLLE style webbing.

Specifications:

- Case Closed: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 4 lb 6.5 oz

SQUAD KIT (CCRK)

COLOR	ITEM#	NSN#
BLK	80-0037	6545-01-543-5721
ODG	80-0038	6545-01-543-3159
COY	80-0039	6545-01-543-2956
MTC	80-0251	----

R_X ONLY

Available Colors:



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

SQUAD KIT

REV081623

Features:

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- Multiple-load carriage configuration capabilities: fanny pack, over the shoulder, attachments for 3-day assault packs and MOLLE style webbing
- Enhanced Hip Stabilization System for greater security during transport
- Trauma Shear Storage on Front for Quick Access

Kit Contents:

- 1x Bag (CCRK® Squad)
- 5 pr Black Talon® Nitrile Trauma Gloves
- 2x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1x Petrolatum Gauze (3 in. x 18 in.)
- 1x HyFin® Chest Seal
- 2x ARS® Needle Decompression Kit (14G x 3.25 in.)
- 2x C-A-T® (Combat Application Tourniquet®)
- 4x ETD® 6 in. Emergency Trauma Dressing
- 4x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1x ETD® Abdominal Emergency Trauma Dressing
- 2x SAM® Splint II
- 1x NAR® Trauma Shears (7.25 in.)
- 4x Polycarbonate Eye Shield
- 1x Surgical Tape (2 in.)
- 1x Combat Casualty Reference Card
- 2x Triage Card (NAR T2 Tag®)



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**