

WALK® WARRIOR AID & LITTER KIT

REV103123

Official Vehicle Life-Saving Kit for the Department of Defense

Allows first-aid treatment and immediate evacuation of one patient in a non-traditional vehicle, plus it allows multiple other trauma patients to be treated until additional transport arrives

From numerous observations and after-action reviews compiled from Operations Enduring and Iraqi Freedom, we learned that a new focus was needed for vehicle first-aid kits. From Warfighters and Tactical Healthcare Professionals, we learned the focus ought to be decreasing preventable combat death at the point of wounding.

Our engineers took the state-of-the-art Talon® 90C (a quad-fold rigid litter with retractable handle system) and designed a package that holds not only the litter, but also a streamlined pocket-system containing self- and buddy-aid medical supplies. The WALK® gives true definition to the term aid and litter kit. The WALK® also contains our Hypothermia Prevention & Management Kit (HPMK™) to help stabilize patients during transport. Light and compact, the WALK® reduces cubic weight and, thereby, evacuation time.

The WALK® interfaces with any standard or non-standard casualty evacuation asset and is ideal as part of the pre-positioning of medical supplies or cache for disaster response. The WALK® has adjustable shoulder straps for occasions when the kit is used to provide treatment and evacuation on target.



Combat Casualty Response Kits®
CCRK®



NORTH AMERICAN RESCUE®


www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**



For the treatment of multiple casualties and the immediate evacuation of one patient

WARRIOR AID & LITTER KIT CONTENTS:



Kit Contents:

- 1x Bag (CCCRK®-WALK®)
- 5 pr Black Talon® Nitrile Trauma Gloves
- 2x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1x Casualty Equipment Bag
- 2x HyFin® Chest Seal
- 2x ARS® Needle Decompression Kit (14G x 3.25 in.)
- 2x C-A-T® (Combat Application Tourniquet®)
- 6x Trauma Emergency Bandage, 6 in.
- 4x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1x ETD® Abdominal Emergency Trauma Dressing
- 2x SAM® Splint II
- 1x NAR® Trauma Shears (7.25 in.)
- 1x Surgical Tape (2 in.)
- 1x GeckoGrip® Multipurpose Tape
- 6x Polycarbonate Eye Shield
- 1x Combat Casualty Reference Card
- 2x Combat Casualty Card (Triage)
- 1x Aviation Panel (Recognition, Orange, 35 in. x 29 in.)
- 1x Talon® II 90C Collapsible Handle Litter
- 1x Hypothermia Prevention & Management Kit™ (HPMK®)
- 4x Tie Down Strap (Universal Litter)

WALK® ReSupply Kit
(Item# 85-0020):
Contains WALK
components for resupply
(excluding nylon bag & litter)



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

WALK® WARRIOR AID & LITTER KIT

Combat Casualty Response Kits®



REV103123



Includes Talon® II 90C quad-fold litter

HPMK-I sold separately

Features:

- Official US Army vehicle life-saving kit
- For treatment of multiple trauma patients and immediate evacuation of one patient
- Includes Talon® II 90C quad-fold litter
- Hypothermia Prevention & Management Kit™ helps stabilize patients during transport
- Interfaces with any standard or non-standard TACEVAC asset
- Accepts both: Talon® II 81C (Retired) & Talon® II 90C Collapsible Handle

Specifications:

- Case Closed: H 22.75 in. x W 16.5 in. x D 11.5 in.
- Weight: 29 lb 15 oz

Available Colors:



| COLOR | ITEM# | NSN# |
|-----------|---------|------------------|
| BLK | 80-0097 | |
| ODG | 80-0060 | |
| COY | 80-0089 | 6545-01-587-1199 |
| RGR | 80-1564 | |
| MTC | 80-1563 | |
| ReSupply* | 85-0020 | |

R_X ONLY

*Contains all WALK components excluding nylon bag & litter



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**