

USCG IFAK

REV081823



Available Colors: **BLK** **COY**

Kit Contents:

- 1 x Bag (USCG -IFAK)
- 2 pr Black Talon® Nitrile Gloves, Large
- 2 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 2 x Petrolatum Gauze (3 in. x 18 in.)
- 1 x HyFin® Chest Seal
- 1 x C-A-T® - Black (Combat Application Tourniquet®)
- 1 x ETD® 6 in. Emergency Trauma Dressing
- 1 x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 1 x Combat Gauze® Z-Fold Hemostatic
- 1 x Polycarbonate Eye Shield
- 1 x GeckoGrip® Rescue Tape



Designed to meet the specific IFAK requirements of the U.S. Coast Guard

Designed to solve problems, not create them, the trim USCG IFAK from North American Rescue can be tucked into small spaces or carried about with ease until it is required for service. Built to the strict requirements of the United States Coast Guard, this improved first aid kit provides for safe first aid until a tactical healthcare professional arrives.

The USCG IFAK is a reliable solution. Made of the highest quality materials, our modular system can be worn by the operator as a drop-down leg rig, hooked directly to the duty belt, or attached to an assault vest or body armor with MOLLE-style connectors.

Features:

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- GhillieTEXTM IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE-style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™
- Sleek Design with no Front Pocket

Specifications:

- Case Closed: H 8.25 in. x W 7.1 in. x D 5.5 in.
- Weight: 2 lb. 10 oz.

COLOR	ITEM #	NSN #
Black	80-0149	6545-01-572-7248
Coyote	80-0158	6545-01-587-8513
ReSupply	85-0104	

R_X ONLY



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**