

# DOOR PANEL KIT

REV082123

## *Easily Accessible and Quick to Deploy*


The Door Panel Kit was created in order to provide a compact, easy to stow major bleeding control kit that can be kept within reach when operating a vehicle. This low profile configuration allows the kit to fit into most common door panel pockets for easy access in the event of an accident or injury on the highway. Kit contents also offer other minor care solutions for stings, smaller lacerations, abrasions, and burns.

### **Features:**

- MOLLE-style attachment system compatible
- 550 cord pulls
- Slim design to fit in door panel storage bin
- Optimized organization of contents

### **Dimensions:**

- Packaged: H 7 in. x W 7 in. x D 1.75 in.
- Weight: 1 lb. 6 oz.



### **Door Panel Kit**

COLOR	GAUZE	ITEM#	NSN#
BLK	Compressed	80-1166	----
RED	Compressed	80-1612	----
BLK	QuikClot®	80-1167	----
RED	QuikClot®	80-0614	----
BLK	CELOX™	80-1615	----
RED	CELOX™	80-1616	----

### **Available Colors:**



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

# DOOR PANEL KIT

REV082123



REAR VIEW



 QTY 1 BLK or ORG C-A-T™ Tourniquet	 QTY 1 Responder Compressed Gauze	 QTY 1 ETD® 4 in. Responder, Flat	 QTY 1 HyFin® Vent Compact, Twin Pk	 QTY 1 Microshield CPR Mask
---	---	---	---	---

 QTY 1 BLK or ORG C-A-T™ Tourniquet	 QTY 1 Choice of Hemostatic	 QTY 1 ETD® 4 in. Responder, Flat	 QTY 1 HyFin® Vent Compact, Twin Pk	 QTY 1 Microshield CPR Mask
--	---	---	---	---

 QTY 1 BurnTec® 4"x4" Dressing	 QTY 3 MED NAP 82K ANTISEPTIC TOWELETTE	 QTY 3 Triple Antibiotic	 QTY 1 BANDAGES KIT	 QTY 3 Safetec Sting Relief
--	---	--	---	---

 QTY 1 BurnTec® 4"x4" Dressing	 QTY 3 MED NAP 82K ANTISEPTIC TOWELETTE	 QTY 3 Triple Antibiotic	 QTY 1 BANDAGES KIT	 QTY 3 Safetec Sting Relief
---	---	--	---	---

 QTY 3 CHAFING OINTMENT	 QTY 1 GEOCKO GRIP FLAT TAPE, 2"	 PAIR 2 RESPONDER BLUE GLOVES, LG	 QTY 1 RESPONDER SCISSORS, LG	 QTY 1 PERMANENT MARKER, LARGE
---	--	---	---	--

 QTY 3 CHAFING OINTMENT	 QTY 1 GEOCKO GRIP FLAT TAPE, 2"	 PAIR 2 RESPONDER BLUE GLOVES, LG	 QTY 1 RESPONDER SCISSORS, LG	 QTY 1 PERMANENT MARKER, LARGE
--	--	---	---	--

 QTY 1 Patch, 2" x 2" Gray/Black
---

## Basic Kit Contents

 QTY 1 Patch, 2" x 2" Gray/Black
--

## w/ Bleeding Control Dressing OR CELOX™ GAUZE