

# RANGE TRAUMA KIT™

REV082223

## Essential equipment to treat injuries common to Shooting and Training Ranges

The Range Trauma Kit by North American Rescue® was designed to address the most common types of injuries that can occur at shooting ranges, training ranges or other locations where firearms are used.

The kit is packaged in a bright orange bag with reflective trim for high visibility. The bag is easy to transport and store. The large central compartment and front Trauma-shear storage offer a versatile, mobile platform with ready-to-use supplies. The kit's contents provide the necessary tools to save lives and to treat injuries by both medically and non-medically trained personnel.

The Range Trauma Kit provides the necessary equipment to treat penetrating injuries, blast injuries, chest wounds, burns, eye injuries, minor lacerations, sprains or fractures and many other common injuries seen on training and shooting ranges.

### Features:

- Designed to address penetrating and other common injuries that occur at shooting ranges, training ranges and other locations where firearms are used
- Bright orange bag with reflective trim for high visibility
- Constructed of durable 1000 Denier Nylon
- Sand & water resistant zipper with silent pulls
- Multiple-load carriage configuration capabilities: fanny pack, over the shoulder, hand-carried, attachments for MOLLE/PALS style webbing
- Trauma Shear Storage on front for quick access

### Specifications:

- Case Closed: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 5 lb. 10 oz.



ORG

## Range Trauma Kit

\*NATO NSN#

COLOR	HEMOSTATIC	ITEM#	NSN#
Orange	ChitoGauze	80-0213	6545-01-610-5411*
Orange	Combat Gauze	80-0299	



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

# RANGE TRAUMA KIT™

REV082223



## Kit Contents:

- 1 x Bag (EMS Rapid Deployment)
- 2 x C-A-T® Tourniquet - Orange
- 1 x Choice of Hemostatic (ChitoGauze® XR PRO or Combat Gauze®)
- 4 x NAR® S-Rolled Gauze
- 4 x NAR® Emergency Trauma Dressing, 6 in.
- 1 x HyFin® Vent Chest Seal (Twin Pk)
- 5 x Sterile Gauze Pads, 4 x 4 in. (2 PK)
- 4 x Polycarbonate Eye Shield
- 2 x Elastic Wrap, 4 in.
- 2 x Triangular Bandage
- 2 x SAM® Splint II
- 1 x NAR® Surgical Tape, 2 in.
- 50 x Flexible Adhesive Bandages
- 2 x Burntec® Burn Dressing, 5 x 10 in.
- 1 x Eye Wash, 4 oz
- 1 x Splinter Forceps, 4.5 in.
- 1 x Responder Shears, Lg.
- 5 pr Responder Gloves - Blue, Lg.
- 1 x CPR Pocket Resuscitator
- 1 x Kit Contents Card



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**