

NUSTAT HEMOSTATIC DRESSING 4 in. x 4 in.

REV082323

Patented dual-fiber technology controls bleeding through absorption at the injury site, which results in a physical aggregation of blood cells and clotting factors. Compact packaging allows for convenient storage and rapid application. Packaged in a convenient vacuum packed pouch, NuStat is available in 4 in. x 4 in. pads and a Tactical Z-fold Dressing, 2 in. x 36 in.

Specifications:

- Patented dual-fiber technology controls bleeding through absorption at the injury site, which results in a physical aggregation of blood cells and clotting factors
- Compact packaging allows for convenient storage and rapid application
- Contains a Single 4 in. x 4 in. Pad in convenient vacuum packed pouch

Dimensions:

- Pad size: 4 in. x 4 in.
- Weight: 1.6 oz

NuStat® Hemostatic Dressing - 4 in. x 4 in.

ITEM #	NSN #
30-0236	----



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®