

S.T.O.R.M. Operator IFAK

Administer Self-Aid/Buddy-Aid to Injured Casualties at the Point of Wounding...

The S.T.O.R.M.® Operator IFAK is the official kit for the Specialized Tactics for Operational Rescue and Medicine (S.T.O.R.M.™) courses developed through a joint effort between the National Tactical Officers Association and the Medical College of Georgia-Center of Operational Medicine.

The S.T.O.R.M.® Operator IFAK is one of our smallest Improved First Aid Kits (IFAK) that was specifically designed to mitigate the three leading causes of preventable death in the tactical environment for penetrating trauma (60% extremity bleeding, 33% tension pneumothorax, 6% airway obstruction).

This kit provides the necessary equipment for the user/operator to administer immediate self-aid/buddy-aid to injured casualties at the point of wounding. Utilizing Tactical Combat Casualty Care (TCCC) Guidelines and lessons learned from the battlefield, this kit provides law enforcement officers and other public safety personnel with the needed tools to reduce the unnecessary loss of life due to the preventable causes of death encountered in the tactical environment. Made to meet the demands of the austere environment, our modular system can be worn as a drop down leg rig, attached directly to the duty belt or attached to the assault vest or body armor with MOLLE/PALS-style connectors.

Features:

- One of NAR's most compact IFAK profiles
- Made of durable CORDURA® Nylon
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE/PALS-style Attachment System Compatible
- Removable Belt and Leg Straps
- Leg Strap with Quick-Grip™

Available Colors:



REV082523



S.T.O.R.M. Operator IFAK

COLOR	ITEM #	NSN #
Black	80-0133	
Coyote	80-0135	
DUC	80-0136	
OD Green	80-0134	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277


PRODUCTS

WITH A

MISSION®

S.T.O.R.M. Operator IFAK

REV082523



Kit Contents:

- 1 x Bag (STORM®)
- 2 x Pair Black Talon® Nitrile Trauma Gloves
- 1 x Nasopharyngeal Airway 28F with Lubricant
- 1 x HyFin® Vent Chest Seal Twin Pack
- 1 x ARS® Needle Decompression Kit (14 G x 3.25 in.)
- 1 x C-A-T® (Combat Application Tourniquet®)
- 1 x ETD™ 6 in. Emergency Trauma Dressing
- 2 x S-Rolled Gauze™ (4.5 in. x 4.1 yd)
- 1 x Polycarbonate Eye Shield
- 1 x Combat Casualty Card
- 1 x Permanent Marker

Dimensions:

- Case Closed:
H 7.4 in. x W 6.9 in.
x D 5 in.
- Weight: 1 lb 12 oz



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®