

# LUMBAR FIRST AID KIT (L-FAK)

REV082523

One of our most compact Individual First Aid Kits on the market today, the L-FAK (Lumbar First Aid Kit) is designed for wear in one of the least utilized areas on the body, the small of the back. This kit empowers First Responders with rapid access to life saving medical equipment without taking up MOLLE real-estate on the wearer's front or sides that are needed for priority access items such as tasers, ammo, or communications equipment.

The L-FAK, one of North American Rescue's most compact personal first aid kit platforms, was designed to carry the necessary items needed for the treatment of penetrating, blast, or other traumatic injuries in the field. This kit comes loaded with critical medical equipment requested by First Responders operating in the line of duty. Treatment content options include Basic, Basic with Combat Gauze, and Advanced versions. NOTE: C-A-T® tourniquet sold separately.

The L-FAK's compact, 420 Denier nylon inner carrier easily opens with the pull of a tab for rapid one-handed access to its contents. Multiple elastic loops on the inner carrier allow for secure organization of medical gear. The 500-denier exterior carrier is equipped with MOLLE style attachment for easy carriage.



## Specifications:

- Durable 500-denier exterior carrier
- Compact design to reduce space needed for life-saving equipment
- Multiple interior elastic loops for securing gear
- Rapid access to medical gear

## Dimensions:

- Packaged - H 4.75 in. x W 5.75 in.
- Basic Weight: 0.40 lbs, w/Combat Gauze Weight: 0.88 lbs, Advanced Weight: 0.94 lbs



## Small of Back First Aid Kit

ITEM #	COLOR	OPTION	NSN #
85-3627	Black	Basic	--
85-3633	Coyote	Basic	--
85-3630	Ranger	Basic	--
85-3628	Black	with Combat Gauze	--
85-3634	Coyote	with Combat Gauze	--
85-3631	Ranger	with Combat Gauze	--
85-3629	Black	Advanced	--
85-3635	Coyote	Advanced	--
85-3632	Ranger	Advanced	--



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
 WITH A  
**MISSION®**

# LUMBAR FIRST AID KIT (L-FAK) .....

REV082523

## Kit Contents:

- 1 x NAR Wound Packing Gauze
- 1 x HyFin® Compact Chest Seal - Twin Pack
- 1 x 4 in. ETD™ Mini Dressing
- 2 x Bear Claw Glove Pack - Tan, Large

### Basic Kit also includes:

- 1 x NAR Wound Packing Gauze

### Kit with Combat Gauze also includes:

- 1 x QuickClot® Combat Gauze Z-Fold

### Advanced Kit also includes:

- 1 x QuickClot® Combat Gauze Z-Fold
- 1 x 28F Nasopharyngeal Airway w/ Lubricant
- 1 x 10ga ARS® Needle



Kits Options also Include: .....



## Small of Back First Aid Kit

ITEM #	COLOR	OPTION	NSN #
85-3627	Black	Basic	--
85-3633	Coyote	Basic	--
85-3630	Ranger	Basic	--
85-3628	Black	with Combat Gauze	--
85-3634	Coyote	with Combat Gauze	--
85-3631	Ranger	with Combat Gauze	--
85-3629	Black	Advanced	--
85-3635	Coyote	Advanced	--
85-3632	Ranger	Advanced	--



NORTH AMERICAN RESCUE®

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

PRODUCTS

WITH

MISSION®