


EZ-IO[®] Manual Extremity Needle Set (Box of 6)

REV082823



38.5 mm, 15G Needle Catheter

ONLY AVAILABLE IN BOX OF 6

**Caution: Federal (USA) law restricts this device to sale by or on the order of a practitioner licensed to use or order use of the device or Federal (USA) law prohibits dispensing without a prescription.*

Features:

- Six different FDA cleared sites for easily accessible iZertion (Proximal Tibia, Proximal Humerus, and Distal Tibia)
- Fluids, medications and drugs delivered immediately to central circulation
- Low profile EZ-Connect included with every needle
- Device is silent and well adapted for confined space
- Body armor can remain in place if needed with manual needle
- Requires minimal pressure
- Product comes in a box of six

Allows The Clinician Complete Control to Immediate Vascular Access...

EZ-IO[®] provides the medical professional immediate vascular access to the central circulation within seconds, delivering medications, intravenous fluids and blood products to adult and pediatric patients. With a specially designed cutting IO needle, manual driver, or small power driver, the EZ-IO[®] allows the clinician complete control – avoiding the use of force. EZ-IO[®] provides rapid, smooth entry into the bone's medullary cavity, creating an immediate conduit to the central circulation.

The EZ-IO[®] Intraosseous Infusion System is the leading choice for intraosseous access across multiple health care settings, and is used by 90 percent of US advanced life support ambulances, over half of US Emergency Departments, and the US Military.

The EZ-IO[®] Manual Needle Set allows medical personnel to immediately secure vascular access manually without removing body armor and exposing the operator to additional risk of injury. Six (6) FDA-cleared peripheral sites provide multiple options for the medic, helping to ensure rapid care for patients.

Specifications:

- Packaged: L 10 in. x W 8.5 in.
- Weight: 20.7 oz.

EZ-IO[®] Products are restricted to Federal Purchases only. Please contact Customer Service for details.

EZ-IO[®] Manual Extremity Needle Set

ITEM#	NSN#
30-0078	----

R_X ONLY



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION[®]**