

NAVAL BOAT RESPONSE AID KIT

REV082923

Critical Life Saving Equipment for Emergency Care at Sea



Boat Response Aid Kits from North American Rescue provide first responders with critical life saving equipment for rendering emergency care at sea in the austere conditions that exist on afloat vessels – both military and commercial.

Kits can be secured using an anti-pilferage seal and the water tight integrity helps to keep your medical products protected.

Meets the rigorous requirements of the U. S. Naval Ships, Military Sealift Command Ships and afloat commercial vessels.

Features:

- Small, compact and rugged hard case
- Airtight, Waterproof and Sealable
- Contains 4 pouches for medical supplies; Bleed, ABD/Burn, Splint, and Eye
- Contains 1 utility pouch which is mounted to the lid of the box
- The pouches were designed to easily and quickly access the medical equipment necessary to respond to specific life threatening medical emergencies
- The pouches are attached by hook and loop to allow the first responder to remove the specific pouch to provide on scene care to the patient

Specifications:

- H 6 in. x W 14 in. x D 10 in.
- Weight: 8 lb 10 oz

Naval Boat Response Aid Kit

ITEM #	NSN
80-0420	6545-01-648-3228



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

NAVAL BOAT RESPONSE AID KIT

REV082923

Kit Contents:

- 1x Naval Response Boat Box, Yellow
- 2x C-A-T® (Combat Application Tourniquet®) Black
- 1x NAR® Emergency Trauma Dressing (ETD®), 6 in.
- 1x NAR® Abdominal/Stump ETD®
- 1x CPR Mask
- 1x CPR Mask, Microshield
- 1x Nasopharyngeal Airway, 28F
- 1x Surgical Lube, 5.0 gm
- 1x HyFin® Vent Chest Seal - Twin Pack
- 1x SAM® Splint
- 1x Cravat, Dry Sterile Burn Dressing, Super XL
- 1x NAR® Elastic Bandage, 3 in.
- 1x ACE Cervical Collar
- 1x Blizzard Blanket, OD Green
- 1x Fox Aluminum Eye Shield w/garter
- 1x Gauze Pads - 2 Pack, 4 in. x 4 in.
- 1x NAR® Trauma Shears, 7.25 in.
- 1x Biological Personal Protection Kit

**Perfect for both
military and
commercial
vessels**



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**