

NAR[®] ELASTIC WRAP BANDAGES

REV082923

Available in a range of widths...

Ideal for low compression and breathable support. Rated for general use and support, elastic bandages from North American Rescue[®] are available in a range of widths and do not contain natural rubber latex.

Designed to wrap snugly around an injured limb or area for support while allowing the area to breathe.

Features:

- Latex-free Elastic Wrap Bandage designed to provide support and compression
- Hook Closure enables the bandage to be wrapped easily and secured without clips
- Ideal for securing hot/cold packs, splinting, sprains, strains, swelling and general soreness
- Available in 4in. flat configuration for reduced cube space

Dimensions:

- **2 inch:** (#30-0027) 2 in. W x 5 ft. L (5 yd. stretched)
Weight: 2.5 oz.
- **3 inch:** (#30-0028) 3 in. W x 5 ft. L (5 yd. stretched)
Weight: 2.5 oz.
- **4 inch:** (#30-0026) 4 in. W x 5 ft. L (5 yd. stretched)
Weight: 2.5 oz.
- **4 inch flat:** (#30-0308) 4 in. W x 5 ft. L (5 yd. stretched)
Weight: 2.5 oz.
- **6 inch:** (#30-0007) 6 in. W x 5 ft. L (5 yd. stretched)
Weight: 3 oz.



Elastic Wrap Bandages

SIZE	ITEM #	NSN #
2 INCH	30-0027	6510-00-935-5820
3 INCH	30-0028	6510-01-230-8702
4 INCH	30-0026	6510-01-171-5130
4 INCH FLAT	30-0308	----
6 INCH	30-0007	6510-01-532-6522



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A **MISSION[®]**