

TCCC IFAK SKILLS KIT *All Combatant*

REV083123

TCCC IFAK Skills Set – All Combatant: The TCCC IFAK Skills Kit – All Combatant is engineered for student use during trauma lanes or scenarios. Each IFAK contains the necessary equipment for initial TCCC care in an easy to attach modular carrying bag is compact that is designed to be worn on the assault vest or body armor via MOLLE-style attachment points.

Features:

- The necessary training equipment for initial TCCC care using the IFAK Skills Set Kit- All Combatant
- Easy to attach modular carrying bag is compact and designed to be worn on the assault vest or body armor via MOLLE-style attachment points

TCCC IFAK Skills Kit Contents:

- 1 x Nylon Bag - Coyote (Tactical Operator Response Kit)
- 1 x C-A-T® Tourniquet - Blue Trainer
- 1 x 6 in. ETD (Flat)
- 1 x NAR S-Rolled Gauze
- 1 x Combat Gauze Trainer
- 1 x Pre-lubed 28F Naso
- 1 x HyFin® Trainer Twin Pack
- 1 x ARS® 14 ga x 3.25 in.
- 1 x Pair Bear Claw Nitrile Gloves - Lg
- 1 x Permanent Marker - Black, Lg
- 1 x Combat Casualty Card
- 1 x Roll 2 in. Surgical Tape



Dimensions:

- H 7 in. x W 4.5 in. x D 2.75 in.
- Weight: 2 lb 8 oz

TCCC IFAK Skills Kit *All Combatant* Rx

OPTION	ITEM #	NSN#
Complete Kit	85-1958	
ReSupply*	80-1028	

*ReSupply includes the contents without the nylon bag



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH  **MISSION®**