

TCCC SKILLS SET *All Combatant*

REV083123

Train up to 25 Students at Individual Skill Stations...

TCCC Skills Set – All Combatant: Provides the necessary equipment to perform required tasks when training at individual skill stations. Includes all the necessary equipment to instruct up to twenty-five (25) students in the treatment of the leading causes of potentially preventable death in the tactical environment. Contents are packaged in a self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory. Addresses extremity hemorrhage, wound packing, chest seal application, needle chest decompression, basic airway adjuncts, splinting, hypothermia prevention & management, and eye injuries.

Features:

- Provides the necessary equipment to train up to twenty-five (25) students at individual skill stations in the treatment of the leading causes of potentially preventable death in the tactical environment
- Packaged in a self-contained, wheeled, rigid case with internal dividers to keep the equipment organized and easy to inventory
- Addresses extremity hemorrhage, wound packing, chest seal application, needle chest decompression, basic airway adjuncts, splinting, hypothermia prevention & management, and eye injuries

Dimensions:

- Hard Storage Case:
External – L 32.5 in. x W 21 in. x D 13 in.
Internal – L 29 in. x W 18 in. x D 10 in.
- Weight: 55 lb


TCCC Skills Set – All Combatant

Kit	Item #	NSN#
Module	80-0968	
ReSupply	80-1029	



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
— WITH —
MISSION®



TCCC SKILLS SET *All Combatant*

REV083123



TCCC Skills Set Contents:

- 1 x Hard Storage Case
- 25 x C-A-T® Tourniquet - Blue Trainer
- 25 x 6 in. ETD (Flat)
- 25 x NAR S-Rolled Gauze
- 25 x Combat Gauze Trainer
- 10 x Pre-lubed 28F Naso
- 12 x HyFin® Trainer Twin Pack
- 10 x ARS® 14 ga x 3.25 in.
- 1 x HPMK®
- 25 x Pair Bear Claw Nitrile Gloves - Lg
- 5 x Permanent Marker - Black, Lg
- 25 x Combat Casualty Card
- 6 x Roll 2 in. Surgical Tape
- 6 x Triangular Bandage - ODG
- 6 x SAM Splint
- 10 x PES, Polycarbonate Eye Shield
- 2 x Sharps Container

#80-1029 ReSupply Kit with contents only is also available



Hard Case Specifications:

- Retractable extension handle
- In-line skate style wheels for easy transport
- Padlock holes include stainless steel padlock protector clips
- Cushion grip handles on three sides for comfortable lifting and carrying
- Waterproof and dust tight design
- Ambient pressure equalization valve
- Molded in hinge latch system
- Trigger release latch system
- Resistant to UV, solvents, corrosion and fungus
- Tan color



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®