

# MARITIME ASSAULT *(BAG ONLY)*

REV090523

## **Sand and Water Resistant Zippers with Silent Pulls**

Split Tail connector and quick grip leg strap attachment or PALS/MOLLE style webbing. This rugged bag can be attached to the operator in three ways:

1. Using the split-tale connector with height adjustment and Quick-Grip™ leg strap
2. PALS/MOLLE style webbing for attachment to the vest/body armor
3. Fasten to the web/gun belt using the quick weave method.

Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration.

## **Specifications:**

- Rapid access and organization is enhanced by elastic storage loops that keep items in place until needed
- 500D CORDURA® Solution-Dyed Nylon with built in Near Infrared Signature Reduction
- Sand and water resistant zippers with silent pulls

## **Dimensions:**

- H 7.4 in. x W 6.9 in. x D 5 in.
- Weight: 8.5 oz

## **Available Colors:**



*(Bag Only)*

## **Maritime Assault Bag (CCRK)**

COLOR	ITEM#	NSN#
BLK	80-0104	8645-01-701-7254
COY	80-0079	----
ODG	80-0105	----



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**