

EXPEDITIONARY JUNIOR MEDIC KIT *EJMK*...

REVI01223

Pouches are designed for quick access to medical equipment necessary to respond to specific life-threatening medical emergencies. Pouches attach via hook-and-loop to allow the first responder to remove the specific pouch needed to provide on-scene medical care to the patient.

Slim-line, low-profile bag is securable using anti-pilferage seals to protect all medical components inside the bag. Kit can be carried as a backpack using shoulder straps that stow away until needed.

Specifications:

- Meets the requirements for the Expeditionary Forces Junior Medic Medical Bag
- Ergo dynamic, low profile & durable casualty response bag constructed of Berry Amendment 500D CORDURA® solution dyed nylon
- GhillieTEX™ hardware for near IR signature reduction and sand/water resistant, self-repairing zippers with silent 550 cord pulls
- Designed to be hung in aircraft or other transport vehicle by top bag handle
- Featuring 6 hook-and-loop clear vinyl internal pouches in various sizes which can be self-labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products

Dimensions:

- H 24 in. x W 15 in. x D 7 in.
- Weight: 20 lb. 3 oz.

EJMK

ITEM #	NSN #
85-2532	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277


PRODUCTS
WITH A **MISSION®**

EXPEDITIONARY JUNIOR MEDIC KIT *EJMK*

REVI01223

EJMK Contents:

- 2 x C-A-T® - Black (Combat Application Tourniquet®)
- 2 x NAR® Compressed Gauze
- 2 x Combat Gauze® Z-Fold Hemostatic Dressing
- 2 x NAR® Emergency Trauma Dressing 6 in.
- 2 x NAR® Emergency Trauma Dressing 4 in.
- 1 x NAR® ETD® ABD/Stump, 6 in.
- 2 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 2 x HyFin® Vent Twin Pack Chest Seal
- 2 x 14G ARS® for Needle Decompression
- 1 x Cyclone® New & Improved BVM
- 2 x 500cc Bags of IV Solution
- 2 x Baxter IV Tubing
- 2 x Survival Blanket
- 2 x SAM® Splint II
- 2 x Triangular Bandage
- 2 x Elastic Bandage, 6 in.
- 2 x Burntec® Burn Dressing, 4 in. x 4 in.
- 2 x Burntec Burn Dressing, 4 in. x 16 in.
- 10 x Antiseptic Towelette
- 10 x Flexible Fabric Bandages
- 2 x Surgical Tape, 1 in.
- 2 x NAR® PPE Kit™
- 3 x Polycarbonate Eye Shield
- 5 pr Bear Claw® Gloves - Sand, Large
- 1 x TORQ tactical light
- 1 x NAR® Sharps Shuttle
- 2 x NAR® Mini Duct Tape, 2 in.
- 1 x NAR® Trauma Shears, Lg
- 1 x Model 7 Rescue Hook
- 1 x Casualty Reference Card
- 2 x Combat Casualty Response Card
- 2 x Permanent Marker - Black, Lg
- 1 x Bag (CCRK® - NAR-5 SAR Aid)



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**