

EODMK Explosive Ordnance Disposal Units Medical Kit

REV091423



The standard casualty response bag meets the requirements for the Explosive Ordnance Disposal Units Medical Bag. Designed to meet the rigorous requirements of the U. S. Air Force Explosive Ordnance Disposal Teams.

Ergo dynamic, low profile & durable casualty response bag constructed of Berry Amendment 500D CORDURA® solution dyed nylon with GhillieTEX® hardware for near IR signature reduction and sand/water resistant, self repairing zippers with silent 550 cord pulls and is designed to be hung in aircraft or other transport vehicle by top bag handle. Featuring 6 hook-and-loop clear vinyl internal pouches in various sizes which can be self labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products.

Pouches are designed for quick access to medical equipment necessary to respond to specific life threatening medical emergencies. Pouches attach via hook-and-loop to allow the first responder to remove the specific pouch needed to provide on-scene medical care to the patient. Slim-line, low-profile bag is securable using anti-pilferage seals to protect all medical components inside the bag. Kit can be carried as a backpack using shoulder straps that stow away until needed.

Features:

- Meets the requirements for the Explosive Ordnance Disposal Units Medical Kit
- Ergo dynamic, low profile & durable casualty response bag constructed of Berry Amendment 500D CORDURA® solution dyed nylon
- GhillieTEX® hardware for near IR signature reduction and sand/water resistant, self repairing zippers with silent 550 cord pulls
- Designed to be hung in aircraft or other transport vehicle by top bag handle
- Featuring 6 hook-and-loop clear vinyl internal pouches in various sizes which can be self labeled for the specific mission set; 4 smaller pouches for hemorrhage control products and 2 larger ones to accommodate airway & intravenous products

Nylon Color: **COY**

Specifications:

- L 24 in. x W 15 in. x D 7 in.
- Weight: 20.3 lbs

COLOR	ITEM#	NSN#
COY	85-3158	




NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

EODMK Explosive Ordnance Disposal Units Medical Kit

REV091423



Kit Contents:

- 1 x Bag (Expeditionary Junior Medical Bag)
- 6 x 14G ARS[®] for Needle Decompression
- 4 x NAR[®] Pre-Lubricated Nasopharyngeal Airway, 28F
- 8 x C-A-T[®] - Black (Combat Application Tourniquet[®])
- 6 x Combat Gauze[®] Z-Fold Hemostatic Dressing
- 4 x CELOX[™]-A Hemostatic Applicator
- 8 x NAR[®] Compressed Gauze
- 4 x NAR[®] Emergency Trauma Dressing, 6 in.
- 2 x NAR[®] Emergency Trauma Dressing, ABD/Stump, 6 in.
- 1 x Burntec[®] Burn Dressing, 5 in. x 10 in.
- 1 x HyFin[®] Vent Chest Seal - Twin Pack
- 1 x HyFin Vent Compact Chest Seal - Twin Pack
- 1 x Cyclone[®] New & Improved BVM
- 2 x NAR[®] Mini Duct Tape, 2"
- 2 x Triangular Bandage
- 1 x Elastic Bandage, 6 in.
- 2 x SAM[®] Splint II
- 1 x NAR[®] Trauma Shears, Lg
- 1 x Model 7 Rescue Hook
- 1 x Hypothermia Prevention and Management Kit (HPMK)
- 1 x Ex-Tract Strap
- 6 x Polycarbonate Eye Shield
- 1 x Rapid Expeditionary Evacuation Litter (REEL)
- 8 pr Bear Claw[®] Gloves - Sand, Lg
- 1 x Survival Blanket
- 1 x First light TORQ with batteries
- 1 x Combat Casualty Reference Card
- 4 x Combat Casualty Response Cards
- 1 x Permanent Marker - Black, Lg
- 4 x Ear Plugs
- 1 x Aviation Signal Panel



NORTH AMERICAN RESCUE[®]
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A **MISSION[®]**