

HIGH RISK WARRANT *Casualty Kit*

REV091823

Rapidly extract casualties to initiate immediate care...

High Risk Warrant Service operations account for the majority of SWAT call-outs nationwide. Due to the rapid dynamic nature of these missions Tactical Health Care Professionals must focus on mobility, adaptation and efficiency.

The High Risk Warrant Casualty Kit™ from North American Rescue allows the Tactical Healthcare Professional the ability to rapidly extract casualties and initiate care for the top three causes of death from penetrating trauma in the tactical environment (60% extremity bleeding / 33% tension pneumothorax / 6% airway obstruction).

Offering a large central compartment and multiple attachment points, this kit maximizes storage space while maintaining a compact profile. Each component is designed to be effective under survival stress and potentially reduce some of the physiological pitfalls associated with it.

Features:

- Multiple carriage capabilities (fanny pack, over the shoulder or attachment with MOLLE-style webbing)
- MOLLE Style Attachments
- Trauma Shears Storage on Front for Quick Access
- Extra Storage Sleeve on Back contains the Dragon Handle System to facilitate Immediate High Threat Extraction
- Durable 500D CORDURA® Nylon
- Provides the capability to address the top three causes of preventable death in the tactical environment

High Risk Warrant Casualty Kit

COLOR	ITEM #	NSN #
Black	80-0121	
OD Green	80-0122	

Available Colors:

BLK **ODG**



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

HIGH RISK WARRANT *Casualty Kit*

REV091823

Kit Contents:

- 2 x C-A-T® - Black (Combat Application Tourniquet®)
- 4 x S-Rolled Gauze (4.5 in. x 4.1 yd)
- 4 x NAR Emergency Trauma Dressing, 6 in.
- 2 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1 x Cyclone® New & Improved BVM
- 1 x Tactical Suction Device
- 1 x HyFin® Vent Chest Seal Twin Pack
- 2 x Petrolatum Gauze (3 in. x 18 in.)
- 2 x ARS® for Needle Decompression (14G x 3.25 in.)
- 5 pr Black Talon® Nitrile Gloves - Black, Lg
- 4x Polycarbonate Eye Shield
- 1 x NAR® Trauma Shears (7.25 in.)
- 2 x Triage Card (NAR T2 Tag®)
- 1 x Dragon Handle System Extraction Device
- 1 x Bag (High Risk Warrant)

Dimensions:

- Case Closed: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 4 lb. 12.8 oz.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**