

M-FAK *Mini First Aid Kit*

REV092424

NAR®'s most compact, versatile Mini First Aid Kit

The M-FAK, North American Rescue®'s Mini First Aid Kit, was designed to be the most compact, multi-use IFAK for delivering immediate critical care for penetrating, blast or other traumatic injuries. At half the size of our popular T.O.R.K.® platform, the M-FAK is currently the smallest platform available in our product offering. Despite its reduced size, this kit comes fully loaded with the critical medical equipment most requested by First Responders operating in the line of duty. Treatment pack out options include Basic, Basic with Combat Gauze® and Advanced versions.



The M-FAK's super compact, 500 Denier nylon bag opens from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 in. internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T® tourniquet. Also included are a loop patch on front of the bag for custom user labeling and (3) three MOLLE strips to attach additional gear such as a TQ/Holder (sold separately) if needed.



M-FAK Mini First Aid Kit

COLOR	CONFIGURATION	ITEM #	NSN#
BLK	Basic	80-0494	
COY	Basic	80-0495	
MTC	Basic	80-0496	
ODG	Basic	80-0497	
GRY	Basic	80-0872	
RGR	Basic	80-0996	
BLK	with Combat Gauze LE	80-0498	
COY	with Combat Gauze LE	80-0499	
MTC	with Combat Gauze LE	80-0500	
ODG	with Combat Gauze LE	80-0501	
GRY	with Combat Gauze LE	80-0873	
RGR	with Combat Gauze LE	80-0998	
BLK	Advanced	80-0502	
COY	Advanced	80-0503	
MTC	Advanced	80-0504	
ODG	Advanced	80-0505	
GRY	Advanced	80-0874	
RGR	Advanced	80-0997	

Nylon Colors:



Features:

- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating penetrating, blast or other traumatic injuries in the line of duty
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access
- Multiple elastic loops for secure gear organization
- Vertical mount can be set to open left-to-right or right-to-left based on shooter preference
- Horizontal mount on a belt allows opening directly to your C-A-T®
- Compatible with NAR® Quick Detach Panels (sold separately)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A MISSION®

M-FAK *Mini First Aid Kit*

REV092424

Content options:

BASIC



with COMBAT GAUZE®



ADVANCED



Specifications:

- L 6 in. x W 3 in. x D 3 in.
- Weight: 13 oz



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**