

GOLDEN HOUR™ MOBILE SERIES CONTAINER PELICAN™

REV092023

...from anywhere in the world — to where they need it most...

The Golden Hour™ Mobile Series Container is a light-weight blood transport capable of carrying four liters of whole blood and keeping the stored blood cool for up to 72 hours, regardless of weather conditions. As the winner of the coveted US Army's Greatest Invention Award, the Series S Container helps deliver whole blood or platelets needed to save lives during the critical first hour after injury. This extremely light-weight container will securely and reliably transport blood and platelets to wherever it is needed most. With Golden Hour™ technology users will experience a reduction to distribution costs, payload risks, and environmental impacts.

Utilizing Thermal Isolation Chamber System (TIC™) technology and vacuum insulation panels, the Golden Hour™ Series S Container will keep the payload isolated from exposure to extreme temperatures with a tight, protective thermal seal for a pre-determined temperature range and duration. Safely deliver whole blood and blood supplies from anywhere in the world — to where they need it most. The containers are easy to assemble, reusable parcels that are constructed of reusable non-toxic materials.

Features:

- Ultra-light storage and transport of four units of whole blood
- Maintains temperature control for 72 hours or more
- Dismount for up to 12 hours without external power
- Equip to make blood in far forward areas accessible
- Capacity: 2L

Dimensions:

- W 9.25 in. x H 8 in. x D 8.25 in.
- Weight: 7.80 lbs



Pelican™ Golden Hour Mobile Series

COLOR	ITEM #
Marine	35-0018
Tan Molle	35-0019



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**