

M-FAK *Mini First Aid Kit - Red*

NAR®'s most compact, versatile Mini First Aid Kit

The M-FAK, North American Rescue®'s Mini First Aid Kit, was designed to be the most compact, multi-use IFAK for delivering immediate critical care for penetrating, blast or other traumatic injuries. At half the size of our popular T.O.R.K.® platform, the M-FAK is currently the smallest platform available in our product offering. Despite its reduced size, this kit comes fully loaded with the critical medical equipment most requested by First Responders operating in the line of duty. Treatment pack out options include Basic, Basic with Bleeding Control Dressing and Advanced versions.

The M-FAK's super compact, 500 Denier nylon bag opens from the side in a clam shell configuration utilizing (2) two main sleeves that open on both ends for easy access and features multiple elastic loops for secure gear organization. The bag's exterior includes a MOLLE backing for mounting in the standard vertical position as well as a 3 in. internal sleeve for accommodating horizontal attachment to your duty belt. The vertical mount can be set to open left-to-right or right-to-left based on shooter preference, while the horizontal mount on a belt allows opening directly to your C-A-T® tourniquet. Also included are a loop patch on front of the bag for custom user labeling and (3) three MOLLE strips to attach additional gear such as a TQ/Holder (sold separately) if needed.



Features:

- Contains the First Responder's most requested critical point-of-wounding medical equipment for treating penetrating, blast or other traumatic injuries in the line of duty
- Super compact, rugged nylon platform that allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Clam shell configuration utilizes (2) two main sleeves that open on both ends for easy access
- Multiple elastic loops for secure gear organization
- Vertical mount can be set to open left-to-right or right-to-left based on shooter preference
- Horizontal mount on a belt allows opening directly to your C-A-T®
- Half the size of our popular T.O.R.K.® platform

**Nylon
Color:** RED

Specifications:

- L 6 in. x W 3 in. x D 3 in.
- Weight: 13 oz

M-FAK Mini First Aid Kit - Red

LEVEL	ITEM#	NSN#
Basic	80-0890	
with Bleeding Control Dressing	80-0891	
Advanced	80-0892	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

M-FAK *Mini First Aid Kit - Red*

REV092023

Content options:


BASIC



with QUIKLOT® BLEEDING CONTROL DRESSING



ADVANCED



*Showing Basic w/
Bleeding Control Dressing®*



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**