

# TACTICAL RAPID DEPLOYMENT KIT

REV092023

## ***Allows reliable and versatile interventions wherever the need arises...***

The Tactical Rapid Deployment Kit is designed to provide personnel with the highest quality trauma supplies for the effective treatment of casualties in unconventional environments.

This kit and its contents are compact, durable and allow reliable and versatile interventions wherever the need arises. Offering a large central compartment and multiple attachment points, this kit maximizes storage space while maintaining a compact profile.

The kit can be carried around the waist as a belly-bag or slung over the shoulder in bandolier fashion. The level of careful detail in its design and configuration sustains the science and experience gained from the valuable lessons learned in the streets and on the urban battlefield.

### **Features:**

- Offers mobility during dynamic events
- Incorporates lessons learned from urban battlefields
- Trauma shears storage on front for quick access
- Designed for efficient use under survival stress
- MOLLE-style attachments

### **Dimensions:**

- Case Closed: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 4 lb 12.8 oz

### **Tactical Rapid Deployment Kit**

COLOR	ITEM #	NSN #
Black	80-0201	
Coyote	80-0203	
DUC	80-0204	
OD Green	85-0202	



### **Available Colors:**



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**



# TACTICAL RAPID DEPLOYMENT KIT

REV092023



*Designed to provide personnel with the highest quality trauma supplies for the effective treatment of casualties in unconventional environments*

## Kit Contents:

- 2 x C-A-T® (Combat Application Tourniquet®) - Black
- 4 x NAR® S-Rolled Gauze
- 4 x NAR® Emergency Trauma Dressing, 6 in.
- 1 x NAR® Abdominal Emergency Trauma Dressing
- 2 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1 x Cyclone® New & Improved BVM
- 1 x Tactical Suction Device
- 1 x HyFin® Vent Chest Seal - Twin Pack
- 2 x SAM® Splint II
- 5 pr Black Talon® Nitrile Trauma Gloves - Black, Lg
- 1 x NAR® Trauma Shears, Lg
- 2 4 x Polycarbonate Eye Shield
- 1 x NAR® Surgical Tape, 2in.
- 2 x Triage Card (NAR T2 Tag®)
- 1 x Bag (Squad)



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**