

MINI MEDIC KIT

REV092123



Breaching
Equipment
Sold
Separately

Nylon Colors:



Compact, multi-purpose trauma aid bag for ultimate versatility at the point of wounding...

The Mini Medic Kit is a compact trauma aid bag designed to provide various capabilities for different mission profiles. External pouches and side/bottom MOLLE allow for external tourniquet or other attachment options. Large pouch on back panel includes a drop-down pouch with abrasive resistant padded bottom allows storage of compact breaching tools. External expansion zipper allows for increased load capacity if needed. The inside of the bag has two small and one large, clear tear-away pouches for medical or other equipment supply storage. Along with two mesh inner pouches and an internal pouch on the inner lid, many storage options exist. The Mini Medic Kit provides a compact, versatile bag equipped with the essential medical equipment to treat multiple casualties.

Features:

- Compact, multi-purpose trauma aid bag with vital medical equipment for point of wounding treatment of multiple injuries
- Two external pockets with two external tourniquet holders
- One external pouch with additional drop-down pouch for use with breaching tools or storage for a hydration pack
- Expansion zipper to increase load capacity
- External MOLLE on sides and bottom of bag / Internal MOLLE rows on inside back of bag
- Three clear tear away pouches on the inside main compartment
- Two internal mesh pouches and one larger internal pouch on inner lid
- Internal tape roll strap

Specifications:

- H 14 in. x W 11 in. x D 5 in.
(Depth expandable to 7.25 in.)
- Weight: 12.5 lb.

Mini Medic Kits

COLOR	LEVEL	ITEM#	NSN#
Black	Basic	80-0979	
Coyote	Basic	80-0980	
OD Green	Basic	80-0981	
Ranger Green	Basic	80-1003	
MultiCam	Basic	80-1231	
Black	Combat Gauze	80-0985	
Coyote	Combat Gauze	80-0986	
OD Green	Combat Gauze	80-0987	
Ranger Green	Combat Gauze	80-1005	
MultiCam	Combat Gauze	80-1232	
Black	Advanced	80-0982	
Coyote	Advanced	80-0983	
OD Green	Advanced	80-0984	
Ranger Green	Advanced	80-1004	
MultiCam	Advanced	80-1233	

Rx ONLY



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

MINI MEDIC KIT

REV092123

QTY 2 C-A-T® TQ Holder	QTY 4 C-A-T® Tourniquet	QTY 4 Gauze S-Rolled	QTY 4 ETD® 4 in. Dressing
QTY 2 HyFin® Vent Chest Seal, Twin Pk	QTY 2 BurnTec® 4"x4" Dressing	QTY 2 SAM® Splint II	QTY 1 NAR® QuikLitter LITE
QTY 2 NAR® Flat Tape, 2"	QTY 1 Gecko Grip® Tape, 2" Mini	QTY 1 NAR Trauma Shears, Lg	PAIR 5 Nitrile Gloves, Lg
QTY 1 Scissor Leash™	QTY 2 Permanent Marker, Large	QTY 4 Combat Casualty Card	QTY 1 Casualty Reference Card
		QTY 1 Patch 2" x 2" Medical Cross	



Designed to provide various capabilities for different mission profiles...

PLUS these Individual Components based upon configuration of choice:

BASIC		COMBAT GAUZE		ADVANCED					
	QTY 4 Gauze S-Rolled		QTY 4 Combat Gauze® Z-Fold		QTY 4 Combat Gauze® Z-Fold	QTY 1 BAC-Pack Bougie-Aided Cric	QTY 2 Pre-lubed NAR® Naso, 28F	QTY 1 Cyclone Pocket BVM	QTY 2 14G ARS® Needle Decomp.



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**