

PAX COOLER BAG *for the Crēdo ProMed Series Four 2 Liter Container*

REV100423

Introducing the North American Rescue® PAX Cooler Bag for the Crēdo ProMed Series Four 2 Liter Container. Carry both your temperature controlled medical supply needs as well as your blood & IV fluid warming device in the same bag for rapid access.

The NAR® Pax Cooler Bag features a main compartment custom sized to fit the Thermal Insulation Chamber - TIC™ and VIP components from the Crēdo ProMed Series Four 2 Liter Container. Simply transfer the removable TIC™ and VIP's into the NAR PAX Blood Cooler Bag from your Crēdo.

In addition to the temperature control compartment, the bag features a detachable section that is ideal for storing the Quantum Blood & Fluid Warming System and other essential IV access equipment.

Specifications:

- Features the non-porous PAX Tech material that provides an easy-to-clean and decontaminate exterior
- Has a detachable section that is ideal for storing the Quantum Blood & Fluid Warming System and other essential IV access equipment
- Compatible with the Thermal Insulation Chamber - TIC™ and VIP components from the 2 Liter Crēdo cooler panels

Dimensions:

- H 7.5 in. x W 11.5 in. x D 10 in.
- Weight: 2 lbs 1.6 oz

Pax Cooler Bag for the Crēdo Promed Series Four 2 Liter Container

ITEM #	NSN #
35-0032	

Use the quick release buckles to detach the front pocket from the bag. The NAR PAX Bag is easy to transport using the reinforced nylon handles or the built-in shoulder strap.

The NAR PAX Bag features the non-porous PAX Tech material that provides an easy-to-clean and decontaminate exterior when compared to traditional woven nylon materials.

Be prepared to easily transport Blood and IV fluids for field deployment with the Blood Cooler Bag for the Crēdo ProMed Series Four 2 Liter Container.



Crēdo ProMed™
sold separately



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**