

# NAR® CLS KIT *Combat Lifesaver Kit*

REV110323

## **The smart alternative...**

The NAR® CLS Kit from North American Rescue® is the smart alternative to the US Army Combat Lifesaver (CLS) Kit that sets the bar for mission specific casualty care platforms. By identifying the needs of the end user's skill set, this CLS kit's innovative design is able to capitalize on mobility, versatility and efficiency. Offering a large central compartment and multiple attachment points, the NAR® CLS Kit comes equipped with optimal medical gear to treat 2 to 4 casualties and is maximized for storage space while maintaining a compact profile. Each component is designed to be effective under survival stress and potentially reduce some of the physiological pitfalls associated with it.

## **Features:**

- Equipped with optimal medical gear to treat 2 to 4 casualties
- Constructed of Durable 1000 Denier Cordura® Solution-Dyed Nylon
- Large central compartment and multiple attachment points
- Sand and water resistant zippers with silent pulls
- ITW Nexus Hardware

## **Specifications:**

- Packaged: L 14 in. x W 14 in. x D 5 in.
- Weight: 4 lb. 14 oz.

## **NAR® CLS KIT**

ITEM #	NSN #
85-1764	



**Available  
Color:**

MTC



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

# NAR® CLS KIT *Combat Lifesaver Kit*

REV110323



## Squad Responder Kit Contents:

- 2 x Combat Application Tourniquet, (C-A-T®) - Black
- 2 x Trauma Dressing, 6 in.
- 2 x NAR® Compressed Gauze, 4.5 in. x 4.1 yd.
- 2 x Combat Gauze® Z-Fold Hemostatic
- 1 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 2 x HyFin® Vent Twin Pack
- 2 x 14G ARS® for Needle Decompression
- 1 x Blizzard Survival Blanket
- 1 x SAM® Splint, Universal Aluminum, 36 in. x 4.25 in.
- 1 x Elastic Bandage, 6 in.
- 1 x NAR® Surgical Tape, 2 in.
- 4 x Fox Eye Shield with garter
- 4 pr Bear Claw® Gloves - Sand, Large
- 6 x Alcohol Prep Pad, 1 in. x 2.5 in.
- 1 x NAR® Trauma Shears, Large, 7.25 in.
- 1 x Scissor Leash™
- 1 x Model 8 Rescue Hook
- 1 x TORQ Tactical Flashlight
- 2 x Combat Casualty Card
- 2 x Permanent Marker - Black, Large
- 1 x Squad Responder Bag, MultiCam

Bag only #80-0932  
also available at  
**NARescue.com**



## Dimensions:

Packaged: L 14 in. x W 14 in. x D 5 in.

Weight: 4 lbs 14 oz



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**