

# FIELD BLOOD TRANSFUSION KIT *FBTK - Donor Module*

REV111623

***“The advantage of the direct transfusion of human blood in cases of severe haemorrhage which we encounter in the emergencies of military surgery cannot to our minds, be overestimated.”<sup>1</sup>***

*The above comment, from a paper presented by Lieutenant Colonel A. Primrose of the Canadian Army Medical Corps in 1916, remains as pertinent today as during those earlier days of military transfusion medicine...*

The Field Blood Transfusion Kit (FBTK) from North American Rescue® is a comprehensive kit providing all the needed supplies to allow the user to rapidly perform a field blood transfusion. All of the supplies conform with the Joint Theater Trauma System Clinical Practice Guidelines. Equipment provided also allows for verifying the blood types of the donor and recipient as well as the collection of blood samples from the donor for additional surveillance testing.

<sup>1</sup>Primrose A, Ryerson ES. The Direct Transfusion of Blood: Its Value in Haemorrhage and Shock in the Treatment of the Wounded in War. Br Med J. 1916;2:384–386.

## Specifications:

- Comprehensive kit allowing the user to rapidly perform a field blood transfusion
- All supplies are compatible with the Joint Theater Trauma System Clinical Practice Guidelines
- Also allows collecting of blood samples from the donor for additional surveillance testing

## Dimensions:

- L 9.5 in. x W 14 in.
- Weight: 0.75 lbs.



## Kit Contents:

QTY 1 Vacutainer Holder, 21G	QTY 4 6 ml Purple Top Vacutainer	QTY 2 6 ml Red Top Vacutainer	QTY 1 Constricting IV Band	QTY 1 6.5" Marked Beaded Cord
QTY 2 NAR Alcohol Prep Pads	QTY 1 4" x 5" Clear Film Dressing	QTY 2 4" x 4" Gauze Pad 2 PK	QTY 2 Swabsticks, Antimicrobial	QTY 1 550 Cord 10"
QTY 1 Blood Bag Single, 450 ml	QTY 1 Blood Bag Sticker	QTY 1 4.25" Kelly Hemostat (Disp.)	QTY 1 572-EWB ASBP Card	
QTY 1 Surgical Tape, 1 in.	PAIR 1 Nitrile Gloves, Lg	QTY 1 Permanent Marker, Large	QTY 1 Donor Card	

## Field Blood Transfusion Kit - Donor Module

ITEM #	NSN #
83-0001	--

R<sub>X</sub> ONLY



**NORTH AMERICAN RESCUE®**  
www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**