

ENT POCKET EXAM LIGHT

REV112723

New LED lamp provides a focused beam of bright white light for superior illumination. Lightweight, waterproof and easy-to-operate twist on/off light activation (2 AAA batteries included).

Features:

- Bright white LED light
- Lightweight
- Waterproof
- Twist on/off light activation
- Batteries included (2 AAA)

Dimensions:

- L 5 in. x Diameter .75 in.
- Weight: 1.6 oz.



ENT Pocket Exam Light

ITEM#	DESCRIPTION	NSN#
10-0007	Exam Light	6515-01-460-9605
ZZ-0038	Replacement Bulb	6515-01-460-9005



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**