

RESCUE ASSAULT TETHER *RAT Strap*

RESCUE

Quick and Effective Option for Hands-Free Extractions...

The RAT™ Strap (Rescue Assault Tether) was designed as a force multiplier during tactical operations. Decreasing manpower demands and increasing efficiency of rescuers, the RAT™ Strap enables hands-free casualty movement and threat suppression simultaneously.

Critical to this task is the ability to rapidly detach from casualty depending on dynamic of environment, which is accomplished with the quick release shackle.

Ideal for shield-assisted, open-field and indoor-corridor rescues. The innovative Daisy Chain configuration is designed to offer multiple reinforced attachment points accommodating the height of the rescuer and the terrain of the rescue site.

Features:

- Designed to meet NFPA Standards
- 1 in. Tech Tape construction with a tensile strength of 4200 lb
- Reinforced stitching, triple bar-tacked loops
- Quick Release Shackle for rapid detachment
- Hands free drag capability
- Decreases manpower needs in “zeroed in” locations.
- Patents Pending

Dimensions:

- Packaged: H 1 in. x W 4 in. x D 5.5 in.
- Extended Length: L 62 in.
- Minimum Break Strength: 2702 lb.
- Weight: 12.7 oz.



RAT Strap

ITEM#	NSN#
90-0002	6530-01-632-7872



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**