


# FIBER OPTIC LARYNGOSCOPE BLADES MAC

REV011224



The Disposable Laryngoscope Blade features the HEINE Classic+ blade shape of the world's best-selling integral fiber optic blade. With proper use the disposable XP blade avoids cross-infection from one patient to the next.

## Features:

- Hygienic and safe
- Optimum light transmission by oversize light guide
- Extremely bright: Highest level of illuminance with HEINE 2.5V or 3.5V LED illuminance within laryngoscope insert
- Blade shape optimized to ease the passage of the tube
- Patient-friendly with reduced risk of injury or damage to teeth
- Extreme torsional stiffness, practically unbreakable
- Blades do not "freeze" or stick to the tongue at low temperatures
- Compatible with all handles to the ISO 7376 (green standard)

## Laryngoscope Blades - Mac

OPTION	ITEM#	NSN#
#2	10-0012	----
#3	10-0013	----
#4	10-0014	----



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**