

FIBER OPTIC LARYNGOSCOPE BLADES *MILLER*

REV011224



◀ **Miller #0**
L 3.1 in. x W 1.1 in. x D 0.9 in.
Weight: 0.3 oz



◀ **Miller #1**
L 4.1 in. x W 1.1 in. x D 0.9 in.
Weight: 0.5 oz



◀ **Miller #2**
L 6.2 in. x W 1.2 in. x D 0.9 in.
Weight: 0.8 oz



◀ **Miller #3**
L 7.7 in. x W 1.2 in. x D 0.9 in.
Weight: 1 oz

The Disposable Laryngoscope Blade features the HEINE Classic+ blade shape of the world's best-selling integral fiber optic blade. With proper use the disposable XP blade avoids cross-infection from one patient to the next.

Features:

- Hygienic and safe.
- Optimum light transmission by oversize light guide.
- Extremely bright: Highest level of illuminance with HEINE 2.5V or 3.5V LED illuminance within laryngoscope insert.
- Blade shape optimized to ease the passage of the tube.
- Patient-friendly with reduced risk of injury or damage to teeth.
- Extreme torsional stiffness, practically unbreakable.
- Blades do not "freeze" or stick to the tongue at low temperatures.
- Compatible with all handles to the ISO 7376 (green standard).

Laryngoscope Blades - Miller

OPTION	ITEM#	NSN#
#0	10-0008	----
#1	10-0009	----
#2	10-0010	----
#3	10-0011	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**