

MEDIC/LEG RIG KIT

Combat Casualty Response Kits®

CCRK®

REV021524

Designed for the Tactical Medic

Tactical Healthcare Professionals who want rapid access to critical basic and advanced equipment by their side choose the Combat Casualty Response Kit®—Medic/Leg Rig™.

The container/bag is constructed of 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction, Ghillie-Tex IR signature reduction hardware and silent, self repairing nylon zippers with 550 pull cords.

The kit can be employed in three ways

1. Using the split-tale connector with height adjustment and Quick-Grip™ leg straps
2. ALS/MOLLE style webbing for att. to the aidbag/ backpack
3. Fasten to the web/gun belt using the quick weave method.

Belt and leg attachment straps can be removed and reattached as needed to allow maximum flexibility during load carriage configuration.

Features:

- 500D CORDURA® Solution-Dyed Nylon with Near Infrared Signature Reduction
- Ghillie-Tex™ IR signature reduction hardware
- Sand & water resistant, self-repairing zipper with silent 550 cord pulls
- MOLLE Style Attachment System Compatible
- Trauma Shear Storage on Back for Quick Access

Specifications:

- Case Closed: H 9.5 in. x W 9 in. x D 7.75 in.
- Weight: 3 lb 9 oz



Available Colors:



MEDIC/LEG RIG KIT (CCRK)

COLOR	ITEM#	NSN#
BLK	80-0046	6545-01-543-3185
ODG	80-0047	6545-01-543-3174
COY	80-0048	6545-01-543-3170
MTC	80-0249	----

Rx ONLY



NORTH AMERICAN RESCUE®

















www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

MEDIC/LEG RIG KIT

REV021524

Rapid access & organization of critical equipment for basic to advanced medical field care

 QTY 2 C-A-T® Tourniquet	 QTY 2 ETD® 6 in. Dressing	 QTY 4 Gauze S-Rolled	 QTY 2 Pre-lubed NAR® Naso, 28F	 QTY 1 Petrolatum Gauze
 QTY 1 HyFin® Occlusive, Ind	 QTY 2 14G ARS® Needle Decomp.	 QTY 1 BOA® IV Constricting Band	 QTY 2 NAR Saline Lock Kit	 PAIR 5 Black Talon Nitrile Gloves, LG
 QTY 1 NAR Trauma Shears, Lg	 QTY 1 Scissor Leash™	 QTY 2 T2® Tag Triage Card	 QTY 1 Casualty Reference Card	 QTY 4 PES Eye Shield
 QTY 1 2" NAR Surgical Tape				

Kit Contents:

- 5 pr Black Talon® Nitrile Gloves
- 2 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28F
- 1 x Petrolatum Gauze (3 in. x 18 in.)
- 1 x HyFin® Chest Seal, Occlusive Single
- 2 x ARS® Needle Decompression Kit (14G x 3.25 in.)
- 2 x C-A-T® (Combat Application Tourniquet®)
- 2 x ETD® 6 in. Emergency Trauma Dressing
- 4 x NAR® S-Rolled Gauze (4.5 in. x 4.1 yd)
- 2 x NAR® Saline Lock Kit
- 1 x BOA® IV Constricting Band
- 1 x NAR® Trauma Shears (7.25 in.)
- 4 x Polycarbonate Eye Shield
- 1 x Scissor Leash™
- 1 x Surgical Tape (2 in.)
- 1 x Combat Casualty Reference Card
- 2 x Triage Card (NAR T2 Tag®)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**