

WOUND MANAGEMENT SIMULATOR

REV022924

Portable training aid to help learn wound management techniques to control bleeding

The NAR Wound Management Simulator is a rugged part-task trainer for persons to learn and simulate the essential life-saving skills of tourniquet application, wound packing and the effective pressure bandage application. Each simulation stump has multiple wound patterns that include: severe laceration and penetrating trauma (entrance and exit wounds).

These part-task training simulators provide the ability to enhance basic bleeding control training and can be used in classroom settings and field training scenarios. Made of a proprietary silicone that is non-environmentally reactive, these rugged and durable simulators provide a portable training aid to help learn wound management techniques to control bleeding.

Features:

- Multiple wounding patterns to enhance training objectives
- Simulates severe laceration, penetrating trauma (entrance and exit wounds), along with the ability to apply tourniquets
- Proprietary silicone that is rugged and durable and will not break down due to environmental exposure
- Can be used in the classroom or in field training scenarios

Specifications:

- L 15.25 in. x Diameter 6 in.
- Weight: 10.8 lb



Wound Management Simulator

Tone	Item#	NSN#
Light	93-0069	
Deep	93-0410	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**