

SHRAPNEL WOUND IN A BOX SIMULATOR

REV032124

The Wounds in a Box® series offers cost effective, compact solutions for skills training. Easy to transport, yet rugged enough for classroom or field use.

The Shrapnel in a Box® task trainer simulates a soft-tissue wound, used to instruct students in proper wound packing technique. It features a hard bone (fractured) inside of the wound cavity to elevate skill level. The representative skin and tissue of the Shrapnel in a Box® is constructed of high quality silicone rubber, with layers of coloration that provide a depth and realism to the product. The PATENTED construction with SOFTTECH® and ROCKTECH® materials a realistic feel, and is very durable and UV resistant. Blood tubing is integrated into the construction to easily create simulated blood flow. The task trainer is contained in its own storage case.

This product can be connected to a bleeding system to provide simulated blood flow that responds to pressure when properly packed. The soft material and patented design provide a realistic feel, and an enhanced training experience. Manual Blood Supply Pump is included.



Specifications:

- Easily transportable, yet rugged enough for classroom and field use
- Patented materials give this product a realistic feel, while still being durable and UV resistant
- Simulator comes in its own storage case
- Manual Blood Supply Pump is included

Dimensions:

- L 6.69 in. x W 4.63 in. x D 3.88 in..
- Weight: 2 lbs. 9.5 oz.

Shrapnel Wound in a Box Simulator

ITEM#	NSN#
93-0240	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

SHRAPNEL WOUND IN A BOX SIMULATOR

REV032124

Construction

The Shrapnel Wound in a Box® features a hard bone (fractured) inside of the wound cavity to elevate skill level. The representative skin and tissue of the Shrapnel Wound in a Box® Simulator is constructed of high-quality silicone rubber, with layers of coloration that provide a depth and realism to the product. In addition, the SOFTTECH® material has a realistic feel, and is very durable and UV resistant.

Operation

The task trainer is contained in its own storage case. The Shrapnel Wound in a Box® includes manual bleeding system. Trainer reinforces wound packing technique for this unique mechanism of injury. Manual blood supply pump is included.

Maintenance

The device should be rinsed in clean, warm water and allowed to dry before storage. The simulation may be removed from the storage case for cleaning. The device may be washed with mild detergents like dish soap. Allow the device to dry fully before storage to prevent mold or mildew growth.

Manual Blood Supply Pump included



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**