

# NIBP CUFF SIMULATOR

REV032124

***The Wounds in a Box® series offers cost effective, compact solutions for skills training. Easy to transport, yet rugged enough for classroom or field use.***

The North American Rescue Simulation Non-Invasive Blood Pressure (NARS NIBP) Monitor Simulator is designed to replicate a portable wrist-style NIBP cuff. This simulator allows medical responders to gather vital signs without instructor intervention. Onboard controls allow various modes of operation. A North American Rescue Simulation (NARS) control tablet is required for operation (*not included*). Display includes systolic & diastolic pressures (mmHg) and heart rate.



## Specifications:

- Allows medical responders to gather vital signs without instructor intervention.
- Various modes of operation via onboard controls
- Includes systolic & diastolic pressures (mmHG) and heart rate on display
- Includes USB charging cable
- Does NOT include NARS control tablet

## Dimensions:

- *Box Dimensions:*  
L 11.63 in. x W 9.69 in. W x D 4.81 in.
- *Product Dimensions:*  
L 11 in. x W 3.63 in. x D 1.88 in.
- *Weight:* 3 lb. 7.5 oz.

## NIBP Cuff Simulator

ITEM#	NSN#
93-0264	----



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**