

K-9 TACTICAL FIELD KIT

Your K-9 partner would defend you with his life; don't leave his to chance.

The NAR K-9 Tactical Field Kit is designed to enable K-9 handlers to provide lifesaving treatment to injured canines at the point of injury or when veterinary care is not immediately available.

This kit and its contents are compact, durable and allow reliable and versatile interventions wherever the need arises. The kit provides the necessary tools to address the major traumatic injuries common to K-9's operating in the military and law enforcement environment as well as common K-9 emergencies such as bloat, heat injury, and toxic ingestion.

Offering a large central compartment and multiple attachment points, this kit maximizes storage space while maintaining a compact profile. The kit can be slung over the shoulder in bandolier fashion, hand-carried, or worn around the waist as a belly-bag.

Ideal For:

- Enables the handler to deliver lifesaving treatment to his canine at the point of injury
- Addresses the major traumatic injuries, such as penetrating and blast-type wounds, which are common to K9s operating in the military and law enforcement arenas
- Provides the essential equipment to render care to the injured K-9 when veterinary care is not immediately available

Dimensions:

- Case Closed: H 7.75 in. x W 13 in. x D 9 in.
- Weight: 5 lb 2 oz



REV032624

K-9 Tactical Field Kit

COLOR	OPTION	ITEM #
Black	Basic	80-1436
Coyote	Basic	80-1437
OD Green	Basic	80-1438
Black	Basic w/ Combat Gauze®	80-0300
Coyote	Basic w/ Combat Gauze®	80-0301
OD Green	Basic w/ Combat Gauze®	80-0303
Black	Basic w/ CELOX™ Gauze	80-1439
Coyote	Basic w/ CELOX™ Gauze	80-1440
OD Green	Basic w/ CELOX™ Gauze	80-1441
Black	Adv. w/ Combat Gauze®	80-1442
Coyote	Adv. w/ Combat Gauze®	80-1443
OD Green	Adv. w/ Combat Gauze®	80-1444
Black	Adv. w/ CELOX™ Gauze	80-1445
Coyote	Adv. w/ CELOX™ Gauze	80-1446
OD Green	Adv. w/ CELOX™ Gauze	80-1447



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A **MISSION®**


K-9 TACTICAL FIELD KIT

REV032624

All kits include:



+plus BASIC kits also include:



+plus ADVANCED kits also include:



All Kits Include:

- 3 x NAR® S-Rolled Gauze
- 2 x ETD® (Emergency Trauma Dressing), 4 in.
- 1 x ETD® (Emergency Trauma Dressing)- Abdominal
- 1 x HyFin® Vent Chest Seal Twin Pack
- 2 x Petrolatum Gauze (3 in. x 18 in.)
- 1 x Eye Wash Solution
- 2 x Polycarbonate Eye Shield
- 1 x Splinter Forceps
- 1 x SAM® Splint II (36 in. x 4.25 in.)
- 1 x Hydrogen Peroxide (8 oz)
- 1 x Black Flexible, Self-Adherent Bandaging Tape
- 2 x Veterinary Elastic Adhesive Tape, 2 in.
- 2 x Veterinary Elastic Adhesive Tape, 4 in.
- 1 x K-9 Litter Kit
- 1 x Tactical Black Muzzle
- 1 x NAR® Trauma Shears, Large
- 1 x Digital Veterinary Thermometer
- 5 pr Responder Nitrile Gloves - Black, Large

also included in Basic:

- 1 x NAR® S-Rolled Gauze
- OR -
- 1 x Combat Gauze® Z-Fold Hemostatic Dressing
- OR -
- 1 x CELOX™ Z-Fold Hemostatic Dressing

also included in Advanced with Combat Gauze®:

- 1 x Combat Gauze® Hemostatic Dressing
- OR -
- 1 x CELOX™ Hemostatic Dressing
- AND -
- 2 x NAR® Pre-lubed Nasopharyngeal Airway, 28F
- 2 x ARS® for Needle Decompression (14G x 3.25 in.)



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**