


CCK CASEVAC Conversion Kit

REV032624



[2] NAR®
Spine Boards

[10] Adjustable
Nylon Support
Straps

[20] Omegalite
4 Wiregate
Carabiners

[4] Spring Clip
Tie Down Straps

Quick Re-Configuration Capabilities...

The CASEVAC Conversion Kit (CCK) by North American Rescue® is designed to enable the MRAP (Mine Resistant Ambush Protected) family of vehicles to be utilized for CASEVAC during critical operations. The CCK provides an effective means of preparing the personnel transport vehicle for safe and secure casualty evacuation in less than 10 minutes. Its quick re-configuration capability from personnel transport to CASEVAC allows one casualty immobilized on a long spine board to be safely loaded and secured for evacuation.

Features:

- Two (2) long spine boards; one serves as the loading platform inside of the vehicle, the other is the casualty board
- Ten (10) adjustable nylon support straps used for suspending and securing the loading platform inside of the vehicle
- Twenty (20) non-locking snap links used to attach the adjustable nylon straps to the vehicle anchor points and the loading platform
- Four (4) spring clip tie down straps: used to secure the casualty board to the loading platform

Dimensions:

- Dimensions: Spine Board - L 72 in. x W 16 in. x D 1.75 in.
- Adjustable Nylon Support Straps - W 2.5 in. x L 14 to 34 in.
- Non-locking Snap Carabiners - 36 mm (1.4 in.) gate opening
- Spring Clip Tiedown Straps - W 1 in. x L 12 in.
- Payload: 1100 lb (598.96 kg)

CASEVAC Conversion Kit

ITEM #	NSN#
50-0036	6530-01-600-4545



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**