

# K-9 RUF KIT

REVO41624



## ***Slim profile trauma aid pack designed to carry a minimum of two complete sets of IFAK contents along with K9 medical supplies***

The NAR® K9 RUF Kit is a slim profile trauma aid pack designed for the K9 Handler or Tracking Team member. The K9 RUF carries a minimum of two complete sets of IFAK contents along with K9 medical supplies. Designed with two tear away internal panels that hold standard IFAK contents and two tear away clear accessory pouches for additional supplies. The internal kit also has a shock cord panel designed to hold a compact, K9 litter. The kit can be worn with a single bandolier style sling strap as well as comes with a waist strap that can be used to secure kit while tracking (waist strap removeable based on user's needs). The slim profile of the case along with the amount of contents of life-saving trauma equipment makes it a necessity for your handler, K9 and tracking team members.

### **Features:**

- Slim profile
- Can be worn with single bandolier style sling strap as well as with a detachable waist strap
- Two tear away internal panels designed to hold standard IFAK contents
- Two additional tear away clear accessory pouches for additional equipment
- Kit contains life-saving trauma equipment and K9 specific equipment
- Internal shock cord area to hold K9 Litter
- External MOLLE on sides of bag for external tourniquet pouches
- Large panel on exterior of kit for ID patches

### **Specifications:**

- Dimensions: H 18 in. x W 8 in. x D 2 in.
- Kit Weight: 6.3 lb
- Bag Weight: 2 lb

### ***K-9 RUF Kit***

	Level	Item#
<b>BLK</b>	Basic	80-1027
	w/ Combat Gauze	80-1026
<b>R<sub>X</sub> ONLY</b>	w/ Celox Gauze	80-1421
	Adv. w/ Combat Gauze	80-1422
	Advanced w/ Celox	80-1423



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

# K-9 RUF KIT

Can be worn with single bandolier style sling strap as well as with a detachable waist strap

REVO41624

Every K-9 RUF KIT includes these core items:

 QTY 2 C-A-T® Tourniquet	 QTY 1 HyFin® Vent Compact, Twin Pk	 QTY 2 NAR® Gauze S-Rolled	 QTY 2 ETD® 6 in. Dressing	 QTY 2 Petrolatum Gauze
 QTY 8 NAR Alcohol Prep Pads	 QTY 1 NAR® K-9 Litter	 QTY 1 Tactical Muzzle	 QTY 2 PES Eye Shield	 QTY 1 Splinter Forceps, 4.5"
 QTY 1 Veterinary Thermometer	 QTY 1 Veterinary Tape 2 in.	 PAIR 4 Responder Black Gloves, LG	 QTY 1 NAR Trauma Shears, Lg	
 QTY 1 Emergency Survival Wrap	 QTY 1 Permanent Marker, Large	 QTY 1 Patch, 2" x 2" Gray/Black		

**PLUS+**

**S-ROLLED GAUZE**  
also includes:

**COMBAT GAUZE**  
also includes:

**CELOX GAUZE**  
also includes:

**PLUS+ K-9 RUF KIT ADVANCED**

**COMBAT GAUZE**  
also includes:

**CELOX GAUZE**  
AND  
also includes:

 QTY 2 NAR® Gauze S-Rolled	 QTY 2 Combat Gauze® Z-Fold	 QTY 2 Celox™ Gauze Z-Fold	- OR -	 QTY 2 Combat Gauze® Z-Fold	 QTY 2 Celox™ Gauze Z-Fold	 QTY 2 Pre-lubed NAR® Naso, 28F	 QTY 2 14G ARS® Needle Decomp.
--	--	--	--------	---	--	---	--



**NORTH AMERICAN RESCUE®**

www.NARescue.com • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**