

BURNTEC® EMERGENCY BURN KIT

REVO41724

First Aid Burn Dressing for 1st, 2nd, & 3rd Degree Thermal, Chemical, and Electrical Burns

The BurnTec® first aid dressing is a hydrogel dressing that is different from other burn first aid dressings due to its stable-sheeting hydrogel structure. The dressing has been strengthened with a medical-quality non-woven fabric which makes it easy to handle. This stable structure eliminates the risk of spilling or damaging the cooling hydrogel. In other words, the hydrogel stays over the burn area and prevents the gel from spilling off the body surface area such as other creams or gels.

The BurnTec® first aid dressing cools the injured area and ensures a soothing effect for up to 24 hours. The sterile non-stick dressing is absorbent, trapping wound secretion and bacteria within the gel. The dressing is permeable to water vapor and oxygen but impermeable to bacteria which helps protect the wound from external contamination while allowing it to breathe. The material is biocompatible and does not cause allergic reactions. BurnTec® is a modern hydrogel dressing used for a variety of skin injuries ranging from first, second, and third-degree burns (thermal, chemical, & electric). This dressing can also be utilized to provide relief from other skin injuries such as sunburns, abrasions, bruises, swelling, or insect bites. Immediate application of the dressing reduces the risk of hypertrophic scars. The non-stick properties allows for painless replacement of the dressing and it does not destroy granulation or new epithelium.



Specifications:

- Easy to handle without spilling or damaging the hydrogel
- Cool, soothing effect for up to 24 hours
- Impermeable to bacteria
- Biocompatible and hypoallergenic
- Reduces risk of hypertrophic scars
- Award winning technology
- CE marked

Small Kit Dimensions:

- Kit: L 8 in. x W 14 in. x D 3 in.
- Weight: 3 lbs 3.2 oz

Large Kit Dimensions:

- Kit: L 10 in. x W 20 in. x D 3 in.
- Weight: 4 lbs 1.6 oz

BurnTec® Emergency Burn Kit

ITEM #	OPTION
80-1661	Small
80-1662	Large



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A
MISSION®

BURNTEC® EMERGENCY BURN KIT

REV041724

Small Kit Contents



For use in the following situations:

- First, second and third-degree burns (thermal, chemical, electric)
- Sunburns
- Abrasions and skin irritations
- Bruises and swelling (insect bites)
- Other skin injuries

Small Kit Contents:

- 1 x BurnTec® Carrying Case
- 2 x BurnTec® Burn Dressing - 4 in. x 4 in.
- 2 x BurnTec® Burn Dressing - 4 in. x 16 in.
- 1 x BurnTec® Burn Face Mask - 12 in. x 16 in.
- 2 x Kerlix™ Gauze - 4.5 in. x 4.1 yd.
- 1 x NAR® Flat Tape - 2 in.
- 1 x Responder Shears, Large

Large Kit Contents:

- 1 x BurnTec® Carrying Case
- 3 x BurnTec® Burn Dressing - 4 in. x 4 in.
- 1 x BurnTec® Burn Dressing - 4 in. x 16 in.
- 1 x BurnTec® Burn Dressing - 8 in. x 16 in.
- 1 x BurnTec® Burn Dressing - 16 in. x 24 in.
- 1 x BurnTec® Burn Face Mask - 12 in. x 16 in.
- 2 x Kerlix™ Gauze - 4.5 in. x 4.1 yd.
- 1 x NAR® Flat Tape - 2 in.
- 1 x Responder Shears, Large

Large Kit Contents



BurnTec® Emergency Burn Kit

ITEM #	OPTION
80-1661	Small
80-1662	Large



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
 WITH A
MISSION®