

# USMC COMBAT LIFESAVER KIT (CLS Kit)

REV032524

## The Official Combat Lifesaver Kit of the Marine Corps

The USMC CLS Kit from North American Rescue® has set the bar for mission specific casualty care platforms. By identifying the needs of the end user's skill set, the CLS Kit's innovative design is able to capitalize on mobility, versatility and efficiency. Offering a large central compartment and multiple attachment points, the kit maximizes storage space while maintaining a compact profile. Each component is designed to be effective under survival stress and potentially reduce some of the physiological pitfalls associated with it.

### Features:

- Official Combat Lifesaver Kit of the Marine Corps
- Constructed of Durable 1000 Denier Cordura® Solution-Dyed Nylon
- Large central compartment and multiple attachment points
- Sand and water resistant zippers with silent pulls
- ITW Nexus Hardware

### Specifications:

- Packaged:  
L 12 in. x W 10 in. x D 8 in.
- Open Length: 37.5 in.
- Weight: 6 lb 14 oz

### USMC CLS KIT

ITEM #	NSN #
80-0159	6545-01-571-4470

Available  
Color: **COY**



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

# USMC COMBAT LIFESAVER KIT (CLS Kit)


REV032524



## USMC CLS KIT:



Bag only #ZZ-0252  
also available at [NARescue.com](http://NARescue.com)



**USMC CLS ReSupply (#80-0424)**  
Contains items for resupply in a clear resealable plastic bag - excludes nylon bag



## USMC CLS Kit Contents:

- 3 x Combat Application Tourniquet® - Black
- 8 x NAR® Compressed Gauze - 4.5 in. x 4.1 yd
- 3 x QuikClot® z-Fold Combat Gauze® - XL
- 3 x Emergency Trauma Dressing - 8 in. x 10 in.
- 1 x Blast Bandage Dressing
- 1 x NAR® Nasopharyngeal Airway - 28F (117 mm) - Pre-lubricated
- 2 x USMC HyFin® Chest Seal Combo Pk (1 x HyFin® Vent & 1 x HyFin® Chest Seal)
- 1 x NAR® Needleless Saline Lock Kit
- 4 x ARS® Needle Decompression Device - 14 ga x 3.25 in.
- 1 x SAM® Splint, Universal Aluminum - 36 in. x 4.25 in.
- 2 x BurnTec® Dressing - 4 in. x 16 in. (10cm x 40cm)
- 3 x Dry Sterile Burn Cravat Dressing - Super XL
- 2 x Fox Eye Shield with Garter
- 1 x Trauma Shears, Large - 7.25 in.
- 8 x Alcohol Prep Pad - 1 in. x 2.5 in.
- 1 x Model 8 Rescue Hook
- 3 pr Bear Claw™ Gloves - Large
- 1 x Gecko Grip® Multi-purpose Flat Tape (6 pack) - 2 in. x 15 in.
- 5 x Combat Casualty Response Card



**NORTH AMERICAN RESCUE®**  
[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**