

# AED PLUS AUTOMATIC DEFIBRILLATOR

ZOLL

REV071624



## **More than CPR. This Is Real CPR Help®.**

Only half of sudden cardiac arrest victims will initially need a shock, but nearly all require high-quality CPR. ZOLL® offers the only AEDs that support rescuers with integrated text, voice, and visual prompts — Real CPR Help®. Its real-time feedback on compression rate and depth helps rescuers administer high-quality CPR. Please note that this product **includes** medical prescription / authorization with purchase.

## **Smart Investment with a Low Total Cost of Ownership**

Longer-lasting batteries and pads require fewer replacement parts, adding up to considerable cost savings over time and rescue readiness in a cardiac arrest emergency.

## **Intuitive One-piece Electrode Pad Design**

Speed is of the essence when treating a victim of sudden cardiac arrest. ZOLL's CPR-D-padz® electrodes feature a unique one-piece, pull-strip design to ensure fast and accurate placement and to minimize the potential of pads being dropped or stuck together.

## **Designed to be Durable**

AEDs are often subjected to rough handling and extreme conditions during an emergency. The ZOLL AED Plus® is highly resistant to dust and moisture and has an ingress protection (IP) rating of 55, among the highest in the industry. It also meets IEC 68-2-27; 100G for shock and military (MIL Standard 810F).

## **Options to Meet Your Preferred Shock Protocol**

The ZOLL AED Plus is available in both semi-automatic (shock button manually pressed by rescuer if a shockable rhythm is detected) or fully automatic (AED automatically administers a shock if a shockable rhythm is detected).



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

# AED PLUS AUTOMATIC DEFIBRILLATOR

ZOLL

REV071624



## ZOLL AED Plus

ITEM#	NSN#
75-0019	----

### Features:

- **INCLUDES MEDICAL PRESCRIPTION / AUTHORIZATION WITH PURCHASE**
- Provides detailed CPR guidance that supports rescuers with integrated text, voice, and visual prompts
- Longer-lasting 5-year battery and electrode pads reduce AED device maintenance requirements
- Includes ZOLL's CPR-D-padz<sup>®</sup> electrodes, which feature a unique one-piece, pull-strip design to ensure fast and accurate placement
- Designed to be durable during emergencies
- Fully Automatic AED

### Package Contents

- 1 x ZOLL<sup>®</sup> AED Plus Automatic Defibrillator
- 1 x ZOLL<sup>®</sup> CPR-D-padz<sup>®</sup> Electrodes
- 1 x Type 123A Photo Flash lithium manganese dioxide batteries (10 pack)
- 1 x Carrying Case - Black

### Specifications:

- **Dimensions:** L 5.25 in. x W 9.5 in. x D 11.5 in.
- **Weight:** 6 lbs. 11.2 oz.
- **Power:** User-replaceable batteries. 10 Type 123A Photo Flash lithium manganese dioxide batteries
- **Device Classification:** Class III and internally powered per EN60601-1
- **Design Standards:** Meets applicable requirements of UL 2601, AAMI DF-39, IEC 601-2-4, EN60601-1, IEC60601-1-2

*For more detailed information regarding this product (such as design standards or battery life), please refer to the ZOLL AED Plus Automatic Defibrillator Spec Sheet.*



**NORTH AMERICAN RESCUE<sup>®</sup>**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION<sup>®</sup>**