

AED PLUS AUTOMATIC DEFIBRILLATOR

ZOLL

REV051624



More than CPR. This Is Real CPR Help®.

Only half of sudden cardiac arrest victims will initially need a shock, but nearly all require high-quality CPR. ZOLL® offers the only AEDs that support rescuers with integrated text, voice, and visual prompts — Real CPR Help®. Its real-time feedback on compression rate and depth helps rescuers administer high-quality CPR.

Smart Investment with a Low Total Cost of Ownership

Longer-lasting batteries and pads require fewer replacement parts, adding up to considerable cost savings over time and rescue readiness in a cardiac arrest emergency.

Intuitive One-piece Electrode Pad Design

Speed is of the essence when treating a victim of sudden cardiac arrest. ZOLL's CPR-D-padz® electrodes feature a unique one-piece, pull-strip design to ensure fast and accurate placement and to minimize the potential of pads being dropped or stuck together.

Designed to be Durable

AEDs are often subjected to rough handling and extreme conditions during an emergency. The ZOLL AED Plus® is highly resistant to dust and moisture and has an ingress protection (IP) rating of 55, among the highest in the industry. It also meets IEC 608-2-27; 100G for shock and military (MIL Standard 810F).

Options to Meet Your Preferred Shock Protocol

The ZOLL AED Plus is available in both semi-automatic (shock button manually pressed by rescuer if a shockable rhythm is detected) or fully automatic (AED automatically administers a shock if a shockable rhythm is detected).



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

AED PLUS AUTOMATIC DEFIBRILLATOR

ZOLL

REV051624



ZOLL AED Plus

ITEM#	NSN#
75-0019	----

Features:

- Provides detailed CPR guidance that supports rescuers with integrated text, voice, and visual prompts
- Longer-lasting 5-year battery and electrode pads reduce AED device maintenance requirements
- Includes ZOLL's CPR-D-padz[®] electrodes, which feature a unique one-piece, pull-strip design to ensure fast and accurate placement
- Designed to be durable during emergencies
- Fully Automatic AED

Package Contents

- 1 x ZOLL[®] AED Plus Automatic Defibrillator
- 1 x ZOLL[®] CPR-D-padz[®] Electrodes
- 1 x Type 123A Photo Flash lithium manganese dioxide batteries (10 pack)
- 1 x Carrying Case - Black

Specifications:

- **Dimensions:** L 5.25 in. x W 9.5 in. x D 11.5 in.
- **Weight:** 6 lbs. 11.2 oz.
- **Power:** User-replaceable batteries. 10 Type 123A Photo Flash lithium manganese dioxide batteries
- **Device Classification:** Class III and internally powered per EN60601-1
- **Design Standards:** Meets applicable requirements of UL 2601, AAMI DF-39, IEC 601-2-4, EN60601-1, IEC60601-1-2

For more detailed information regarding this product (such as design standards or battery life), please refer to the ZOLL AED Plus Automatic Defibrillator Spec Sheet.



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION[®]**