

AED 3 AUTOMATIC DEFIBRILLATOR

ZOLL

REV071624

The ZOLL AED 3[®] automated external defibrillator (AED) provides detailed CPR guidance with proprietary Real CPR Help[®] technology. More powerful full-color images and a CPR cycle timer, combined with voice and text prompts, help guide responders through the rescue. The Real CPR Help integrated, real-time CPR feedback tells and shows rescuers when they need to “Push Harder” or whether they are doing “Good Compressions,” so they know when they are administering high-quality CPR. Please note that product **includes** medical description / authorization with purchase.

RapidShock™ Analysis

The ZOLL AED 3[®] with RapidShock[®] analysis enables the shortest rhythm analysis for more continuous care and more life-saving CPR. RapidShock analysis allows the ZOLL AED 3 to deliver a shock, if needed, in as little as 5 seconds. Minimizing the pre-shock pause and providing more CPR can improve patient outcomes.*

Research shows that minimizing time to shock after the CPR cycle ends may improve survival.¹ The AHA 2015 guidelines state: Shorter interruptions for shock delivery are associated with greater likelihood of shock success, greater likelihood of return of spontaneous circulation, and greater survival to hospital discharge.

Integrated Pediatric Rescue

Pediatric rescue is made easier with ZOLL's unique CPR Uni-padz[®] electrodes and a child-mode setting. Universal pad design provides rescuers a single solution to treat both adult and pediatric victims of SCA by using the same set of electrode pads and simply activating child mode.

Low Total Cost of Ownership

Longer-lasting 5-year battery and electrode pads reduce AED device maintenance requirements. Less maintenance can improve device readiness and has the added benefit of reducing the total cost of ownership over the life of your ZOLL AED 3 device.

Enhanced Real CPR Help[®]



Optional WiFi Connectivity

Easily manage your suite of devices via the cloud. ZOLL AED 3 devices are equipped with optional Wi-Fi connectivity allow seamless communication with AED program management systems through your wireless network. Once enabled, the ZOLL AED 3 will automatically report its readiness status, report a failed self-test, and will send email alerts so you can take the appropriate corrective action and make sure your AEDs are ready in an emergency.

** On the ZOLL AED 3 Automatic, the pre-shock pause time is 5 seconds longer due to warning to stand clear of patient and shock countdown before delivery.*



NORTH AMERICAN RESCUE[®]

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION[®]**

PRODUCTS
WITH A **MISSION[®]**

AED 3 AUTOMATIC DEFIBRILLATOR

ZOLL

REV071624



ZOLL AED 3

ITEM#	NSN#
75-0024	----

Package Contents

- 1 x ZOLL® AED 3 Automatic Defibrillator
- 1 x CPR Uni-Padz™ Electrodes
- 1 x ZOLL® AED 3 Battery Pack

Features:

- **INCLUDES MEDICAL PRESCRIPTION / AUTHORIZATION WITH PURCHASE**
- Provides detailed CPR guidance that supports rescuers with integrated text, voice, and vivid color visual prompts
- RapidShock® Analysis enables short rhythm analysis for continuous care and more life-saving CPR
- Pediatric rescue is made easier with ZOLL's unique CPR Uni-padz® electrodes and a child-mode setting
- Longer-lasting 5-year battery and electrode pads reduce AED device maintenance requirements
- Equipped with optional Wi-Fi connectivity, allowing seamless communication with AED program management systems through your wireless network

Specifications:

- **Dimensions:** L 5 in. x W 9.3 in. x D 9.7 in.
- **Weight:** 5 lbs. 8 oz. (with battery)
- **Power:** Battery
- **WiFi-enabled Devices:** Wireless 802.11 a/b/g/n; Security Protocols, WPA1, WPA 2, WPA Personal, WPA Enterprise; Internal Clock Synchronization: Coordinated Universal Time (UTC) synchronization occurs when communicating with the ZOLL Online server.
- **Rescue Accessory pack Contents:** Non-latex gloves, a CPR barrier mask, razor, paper towel, and towelette.
- **Device Classification:** Class II and internally powered per EN60601-1
- **Design Standards:** Meets applicable requirements of EN 60601-1, IEC 60601-1-11, IEC 60601-2-4

For more detailed information regarding this product (such as design standards or battery life), please refer to the ZOLL AED 3 Automatic Defibrillator Spec Sheet.



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**