

RIG SERIES EAGLE IFAK

REV092524



Low profile IFAK with rapid access to your C-A-T® tourniquet...

The RIG Series Eagle IFAK is a new low profile IFAK designed with a protective hood featuring a quick latch release to keep your life-saving equipment protected and allows rapid access for deployment.

Highlighted design advancements include easy-one hand quick release hood latch and zippers, detachable inner shock cord panel that detaches from the IFAK bag and attachment system that allows both vertical (MOLLE backing) and horizontal (3 in. belt loop) carriage orientations.

Features:

- Tourniquet attached for quick deployment by banding it to inside of front opening flap
- Protective hood with quick latch release helps protect IFAK components and includes sewn in embroidered medical cross identifier
- Allows attachment both vertically (MOLLE backing) and horizontally (3 in. belt loop)
- Quick pull latch and tab for easy one-handed opening
- Detachable (hook & loop backed) shock-cord inner panel adjusts to secure internal components & is easily removable for rapid access)
- Compatible with NAR® Quick Detach Panels (sold separately)


Specifications:

- Bag Dimensions:
H 8 in. x W 4 in. x D 3 in.
Weight: 6 oz
- Kit Dimensions:
H 8 in. x W 4 in. x D 3 in.
Weight: 1.18 lb

RIG Series Eagle IFAK

Color	Level	Item#	NSN#
Black	Basic	80-0509	
Coyote	Basic	80-0510	
Multi-Cam	Basic	80-0511	
Ranger Green	Basic	80-0991	
Black	Combat Gauze	80-0517	
Coyote	Combat Gauze	80-0518	
Multi-Cam	Combat Gauze	80-0519	
Ranger Green	Combat Gauze	80-0993	
Black	Advanced	80-0513	
Coyote	Advanced	80-0514	
Multi-Cam	Advanced	80-0515	
Ranger Green	Advanced	80-0992	

Nylon Colors:



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

RIG SERIES EAGLE IFAK

REV092524



Every RIG SERIES EAGLE IFAK includes these core items:

 QTY 1 C-A-T® Tourniquet	 QTY 1 Gauze, Wound Packing	 QTY 1 ETD® 4 in. Dressing	 QTY 1 HyFin® Vent Chest Seal, Twin Pk	 PAIR 1 Black Talon® Nitrile Gloves, LG
---	---	--	--	---

PLUS these Individual Components based upon your configuration of choice:

BASIC  QTY 1 Gauze, Wound Packing	COMBAT GAUZE  QTY 1 Combat Gauze® LE	ADVANCED  QTY 1 Combat Gauze® LE	 QTY 1 Pre-lubed NAR® Naso, 28F	 QTY 1 14g ARS® Needle Decomp.
--	--	--	---	--



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**