

NAR® NEEDLELESS SALINE LOCK KIT

REV050324



IV Safety Catheter & Needle Free Valve Port Reduce or Eliminate Accidental Needle-Stick Exposures

North American Rescue®'s Needleless Saline Lock Kit allows for quick and effective IV cannulation under the most demanding conditions. Our convenient kit is packaged in a compact, vacuum-sealed, rugged pouch to keep the contents in operational ready condition until needed for use. The NAR® Needleless Saline Lock Kit contains an IV safety catheter and a needle-free valve port to reduce or eliminate accidental needle-stick exposures to tactical healthcare providers.

Features:

- Provides essential components for establishing a saline lock in the most demanding conditions.
- Contains engineered safety devices to reduce or eliminate healthcare provider exposure to accidental needle sticks.
- Includes a clear, adhesive dressing and the NAR® Raptor IV Securing Device.
- Packaged in a compact, vacuum-sealed, rugged pouch.

Specifications:

- Pkg: H 5.5 in. x W 3.25 in. x D 1 in.
- Weight: 1.3 oz

ITEM#	NSN#
30-0046	6515-01-607-1943



 QTY 1 Constricting IV Band	 QTY 1 Raptor IV Securing Device	 QTY 2 NAR Alcohol Prep Pads	 QTY 1 2" x 2" Gauze Pad 2 PK	 QTY 2 2.4" x 2.75" Clear Film Dressing
 QTY 1 Needle Free Valve Port	 QTY 1 Hypodermic Needle, 18G x 1.5"	 QTY 1 IV Safety Catheter, 18G	 QTY 1 Luer-Lock Syringe, 10cc	



NORTH AMERICAN RESCUE®
www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**