

TCCC MAN[®] SIMULATOR

with NARS Instructor Tablet

Rugged, state-of-the-art TCCC MAN[®] training simulator provides enhanced training in Tactical Combat Casualty Care (TCCC) from Point of Injury to transfer of higher care!

- **All-in-one, hyper-realistic trauma simulation solution** enabling instructors to build **multiple realistic traumatic casualty scenarios** with **reactive patients** ideal for use in classroom settings, trauma lanes and field training scenarios
- Specifically **designed** for the **application of skills** and **education** utilizing protocols for **Tactical Combat Casualty Care (TCCC)**



TCCC MAN

ITEM #

93-0958

Features:

- Remote controlled bleeding from multiple sites
- Arterial and venous bleeding
- Quick-exchange Lithium-ion battery
- Large capacity on-board bleeding bag
- Weight 150 lbs
- Industry-leading durability - withstands drops, falls, and drags
- Welded steel frame, steel joints and bones
- Austere environmental use (rain, snow, mud, sand)

TCCC Interventions Supported:

- Massive hemorrhage control (tourniquet application and wound packing)
- Chest seal application
- Pleural needle decompression (multiple sites)
- Airway management (NP/OP airways, supraglottic devices, endotracheal intubation)
- Tibial intraosseous site



NORTH AMERICAN RESCUE[®]

Simulation

866-839-3455 · NARescue.com

REV050224

Discover more hyper-realistic **Training Simulators** available at **NARescue.com** today!

TCCC MAN[®]

Disconnect your student's dependency on instructor "over-the-shoulder" input

Kit Includes:

- 1 x Remote instructor tablet
- 1 x Blood reservoir pressure fill tank
- 1 x Velcro digital pattern uniform
- 2 x 3-wheeled deployment-style bags for safe storage and transport
- 1 x NAR Training IFAK
- 1 x Operator's Manual

Injuries Included:

- GSW chest with anterior entry and mid-axillary exit wound
- Deep packable pelvic wound
- Complete below the knee amputation leg (TQ responsive)
- Unstable pelvis with palpable crepitus

TCCC MAN

ITEM #

93-0958

MILSPEC INSTRUCTOR TABLET *with Scenario Preset Builder Software*

Features:

- Features proprietary communication link that requires NO Wi-Fi, Bluetooth, or LAN connections
- Simultaneous long-range communications with numerous NARS devices
- NO interference risk with telemetry, avionics, or weapon systems



NORTH AMERICAN RESCUE[®]

Simulation

www.NAR-simulation.com

FROM TRAINING TO DEPLOYMENT - NORTH AMERICAN RESCUE DELIVERS



NORTH AMERICAN RESCUE[®]

Education and Training

www.NAReducation.com



NARMD
MEDICAL DEPOT

www.NARmd.com



NORTH AMERICAN RESCUE[®]

www.NARescue.com

REV050224