

CRIC IN A BOX SIMULATOR

REV052324

The Wounds in a Box® series offers cost effective, compact solutions for skills training. Easy to transport, yet rugged enough for classroom or field use.

The Cric in a Box™ is a portable task trainer, perfect for practicing your skills with cricothyrotomy devices.

The Cric in a Box™ incorporates the clinical landmarks of the Hyoid, Thyroid, and Cricoid cartilages. In addition, the tracheal rings provide tactile confirmation of successful placement.

This part task trainer provides a realistic feeling for training to perform a cricothyrotomy. It comes with simulated skin and membranes that allow the trainee to practice surgical incisions required to establish an emergency airway.

The simulated trachea has interior rings to provide tactile feedback when probed. It is constructed of a tough, yet flexible composite that allows the use of a tracheal hook. The Cric in a Box™ comes with 5 non-bleeding neck skins and 25 membranes.



*** Item#: 93-0103 Replacement Cric Membrane 25-Pack available!**

Specifications:

- Easily transportable, yet rugged enough for classroom and field use
- Patented materials give this product a realistic feel, while still being durable and UV resistant
- Simulator comes in its own storage case
- 5 non-bleeding neck skins and 25 membranes are included

Dimensions:

- L 9.13 in. x W 5.88 in. x D 4.56 in.
- Weight: 3 lbs. 3.9 oz.

Cric in a Box™ Simulator

| ITEM# | NSN# |
|---------|------|
| 93-9066 | ---- |



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**