

MEDIC ASSAULT KIT *BASIC AND ADVANCED*

REV062724

Low Profile, Contoured Fit for Trauma and Rescue Aid

The Medic Assault Rescue Kit is our newest compact, low-profile trauma & rescue aid bag. This bag is unique because it is designed with a middle hinge, allowing a contoured fit for room breaching and clearing.

Two main bag compartments are accessible from each end of bag independently. The front of the bag has two large front pockets with full MOLLE on the exterior. Inside, there are four large mesh inner pouches that are secured to the bag interior with hook & loop attachment and provide organization via removable hook & loop labels for customization. The Basic Rescue kit comes with two Bleed pouches, one Splinting pouch, and one Airway pouch. The Advanced Rescue kit comes with two Bleed pouches, one IV pouch, and one Airway pouch.

MOLLE on back and bottom of the exterior of the nylon bag provides additional carry capability. All silent zipper pulls. Removable backpack straps allow attachment directly to body armor if desired (requires additional straps/hardware – not included).


Available in six 500D Cordura colors, including MTC Alpine.

Features:

- Slim Design and Low Empty Weight
- Middle hinge design to allow contouring to the body decreasing profile
- Large mesh internal bags
- Attachable to body armor

Specs:

- H 22 in. x W 11.5 in x D 5 in.
- Weight: 12.2 lbs. (*Basic*); 12.5 lbs. (*Adv.*)



Rx only

COLOR	OPTION	ITEM#	NSN#
BLK	Basic	80-1650	----
COY	Basic	80-1651	----
ODG	Basic	80-1653	----
RGR	Basic	80-1654	----
MTC	Basic	80-1652	----
MTC ALPINE	Basic	80-1649	----
BLK	Advanced	80-1731	----
COY	Advanced	80-1732	----
ODG	Advanced	80-1734	----
RGR	Advanced	80-1735	----
MTC	Advanced	80-1733	----
MTC ALPINE	Advanced	80-1730	----



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

MEDIC ASSAULT KIT

Includes:



PLUS Advanced Kit Also Includes:



CONTENTS:

- 2 x C-A-T® Tourniquet - Black
- 4 x NAR® Wound Packing Gauze, 3 in. x 5 yd.
- 2 x NAR® Emergency Trauma Dressing (ETD®), 6 in.
- 2 x NAR® Flat ETD®, 4 in.
- 1 x Abdominal/Stump ETD®, 12 in. x 12 in.
- 2 x NAR® Pre-Lubricated Nasopharyngeal Airway, 28FR
- 1 x Tactical Suction Device
- 2 x HyFin® Vent Chest Seal Dressing - Twin Pk
- 2 x NAR® Emergency Survival Wrap
- 1 x Ex-Tract Strap
- 2 x SAM® Splint II
- 2 x Triangular Bandage - OD Green
- 2 x Polycarbonate Eye Shield (PES) - Orange
- 1 x NAR® Trauma Shears - Black, Large
- 1 x Scissor Leash™ - Black
- 1 x Aviation Panel Kit - 35 in. x 29 in.
- 2 x Combat Casualty Card
- 2 x BLEED Nametape - Black/Gold
- 1 x AIRWAY Nametape - Black/Gold
- 1 x IV Nametape - Black/Gold
- 1 x Permanent Marker - Black, Large

PLUS (ADV)

- 2 x Quikclot® Z-Fold Combat Gauze, 3 in. x 4 yd.
- 1 x Cyclone® Pocket Bag Valve Mask (BVM)
- 1 x Bougie Aided Cric (BAC) Pack
- 2 x Enhanced ARS™ Decompression Needle, 10 ga
- 2 x NAR® Saline Lock Kit
- 1 x BOA® IV Constricting Band
- 1 x Intraosseous EZ-IO Manual Needle (Box of 6)
- 1 x NAR® Mini Duct Tape, 2 in.
- 1 x Halstead Straight Mosquito Forceps, 5 in.
- 1 x Sharps Shuttle
- 1 x Armadillo Medicine Case - Black

REV062724



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277