

# MEDIC ASSAULT *BAG ONLY*

REV062724

The Medic Assault Rescue Bag is our newest compact, low-profile trauma & rescue aid bag. This bag is unique because it is designed with a middle hinge, allowing a contoured fit for room breaching and clearing.

Two main bag compartments are accessible from each end of bag independently. The front of the bag has two large front pockets with full MOLLE on the exterior. Inside, there are four large mesh inner pouches that are secured to the bag interior with hook & loop attachment and provide organization via removable hook & loop labels for customization.

MOLLE on back and bottom of the exterior of the nylon bag provides additional carry capability. All silent zipper pulls. Removable backpack straps allow attachment directly to body armor if desired (requires additional straps/hardware – not included).

Available in six 500D Cordura colors, including MTC Alpine.




## Features:

- Slim Design and Low Empty Weight
- Middle hinge design to allow contouring to the body decreasing profile
- Large mesh internal bags
- Attachable to body armor

## Bag Specs:

- H 22 in. x W 11.5 in x D 3 in.
- Weight: 2.2 lbs.



COLOR	ITEM#	NSN#
BLK	80-1523	----
COY	80-1524	----
ODG	80-1526	----
RGR	80-1527	----
MTC	80-1525	----
MTC ALPINE	80-1522	----



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**