

MEDIC ASSAULT RESCUE KIT *BASIC AND ADVANCED*

REV062724

Low Profile, Contoured Fit for Trauma and Rescue Aid

The Medic Assault Rescue Kit is our newest compact, low-profile trauma & rescue aid bag. This bag is unique because it is designed with a middle hinge, allowing a contoured fit for room breaching and clearing.

Two main bag compartments are accessible from each end of bag independently. The front of the bag has two large front pockets with full MOLLE on the exterior. Inside, there are four large mesh inner pouches that are secured to the bag interior with hook & loop attachment and provide organization via removable hook & loop labels for customization. The Basic Rescue kit comes with two Bleed pouches, one Splinting pouch, and one Airway pouch.

Features:

- Slim Design and Low Empty Weight
- Middle hinge design to allow contouring to the body decreasing profile
- Large mesh internal bags
- Attachable to body armor

Specs:

- H 22 in. x W 11.5 in x D 5 in.
- Weight:
12.2 lbs. (*Basic*)
12.5 lbs. (*Adv.*)

The Advanced Rescue kit comes with two Bleed pouches, one IV pouch, and one Airway pouch. Additionally, the Rescue versions of the Medic Assault Kits come with carabiners, rope grabs, assorted ropes, and tubular webbing.

MOLLE on back and bottom of the exterior of the nylon bag provides additional carry capability. All silent zipper pulls. Removable backpack straps allow attachment directly to body armor if desired (requires additional straps/hardware – not included).

Available in six 500D Cordura colors, including MTC Alpine.



COLOR	OPTION	ITEM#	NSN#
BLK	Basic	80-1737	----
COY	Basic	80-1738	----
ODG	Basic	80-1740	----
RGR	Basic	80-1741	----
MTC	Basic	80-1739	----
MTC ALPINE	Basic	80-1736	----
R _X ONLY			
BLK	Advanced	80-1743	----
COY	Advanced	80-1744	----
ODG	Advanced	80-1746	----
RGR	Advanced	80-1747	----
MTC	Advanced	80-1745	----
MTC ALPINE	Advanced	80-1742	----



NORTH AMERICAN RESCUE®

































www.NARescue.com • 888.689.6277

PRODUCTS
WITH A **MISSION®**

MEDIC ASSAULT RESCUE KIT


REV062724

Includes:

 QTY 2 C-A-T® Tourniquet	 QTY 4 NAR® Gauze, Wound Packing	 QTY 2 ETD® 6 in. Dressing	 QTY 2 ETD® 4 in. Dressing	 QTY 1 ETD® Abd./Stump	 QTY 2 Pre-lubed NAR® Naso, 28F	 QTY 1 Tactical Suction Device
 QTY 2 HyFin® Vent Chest Seal, Twin Pk	 QTY 2 Emergency Survival Wrap	 QTY 1 Ex-TRACT Strap	 QTY 2 SAM® Splint II	 QTY 2 Triangular Bandage	 QTY 2 PES Eye Shield	 QTY 1 NAR Trauma Shears, Lg
 QTY 1 Scissor Leash™	 QTY 1 Aviation Panel, 35" x 29"	 QTY 2 Combat Casualty Card	 QTY 2 Screw Lock Carabiner	 QTY 2 Oval Auto- Lock Carabiner	 QTY 2 HMS Auto- Lock Carabiner	 QTY 1 Rope Grab /Ascender
 QTY 1 Anti-Panic Descent Device	 QTY 1 Pulley w/ Rope Grab	 QTY 1 HMPE Cord Sling 120cm	 QTY 1 7mm VT Prusik	 QTY 1 50 ft - 8 mm HAVOC Rope	 QTY 1 75 ft - 6.8 mm RIT 900 Rope	 QTY 1 Nylon Tubular Webbing
 QTY 2 Nametape BLEED	 QTY 1 Nametape AIRWAY	 QTY 1 Nametape SPLINTING	 QTY 1 Permanent Marker, Large			

PLUS Advanced Kit Also Includes:

 QTY 2 Combat Gauze® Z-Fold	 QTY 1 Cyclone® Pocket BVM	 QTY 1 BAC-Pack Bougie-Aided Cric	 QTY 2 Enhanced ARS Needle, 10G
 QTY 2 NAR® Saline Lock Kit	 QTY 1 BOA® IV Constricting Band	 QTY 1 Intraosseous EZ-IO Manual Needle	 QTY 1 NAR® Mini Duct Tape, 2 in.
 QTY 1 5" Halstead Straight Mosquito Forceps	 QTY 1 Sharps Shuttle	 QTY 1 Armadillo Med Case	



NORTH AMERICAN RESCUE®

www.NARescue.com • 888.689.6277