

# ZOLL AED PLUS CPR-D-PADZ®

*One Piece Electrode with Real CPR Help®*

REV071624

ZOLL's CPR-D-padz® Electrodes's one-piece design ensures fast and accurate placement on a patient. This product is for use on adult patients and children over the age of 8, or greater than 55 pounds (25 kilograms). Additionally, Real CPR Help technology is built into the electrodes, which help sense the rate and depth of chest compressions.

Every CPR-D-padz® Electrodes are packaged with Gloves, Barrier Mask, Scissors, Razor, a Wet Wipe and a Dry Wipe. CPR-D-padz® Electrodes have a shelf life of five years.

## Features:

- **INCLUDES MEDICAL PRESCRIPTION / AUTHORIZATION WITH PURCHASE**
- Compatible with the ZOLL® AED Plus Automatic Defibrillator
- For use with adult patients and children over the age of 8 years or greater than 55 pounds (25 kilograms)
- Has integrated Real CPR Help® technology
- Has a five-year shelf life

## Dimensions:

- Package: L 8.5 in. x W 8.25 in. x D 1.75 in.
- Packaged Weight: 13 oz.



## ZOLL AED PLUS CPR-D-PADZ®

ITEM#	NSN#
75-0023	----



**NORTH AMERICAN RESCUE®**

[www.NARescue.com](http://www.NARescue.com) • 888.689.6277

**PRODUCTS**  
WITH A **MISSION®**

**PRODUCTS**  
WITH A **MISSION®**